



# LJ Charts in TTI: Interpretation & Trouble shooting

Ankit Mathur  
Consultant Transfusion Medicine &  
Transplant Immunology



Rotary Bangalore  &  Initiative

# Quality Control: Essentials

- Policies & Procedure: Strict adherence
- Training of staff
- Documentation
- Each test run must include one **full set of controls** that yield results within the limits of standard for acceptability and validity
- **Physical parameters** (incubation time, temperature range, reagents conc.) should be strictly followed
- Use test kits **before expiry** date

# Controls

## Internal controls

- **Set of controls (Positive & Negative) provided along with the kit**
- **To be used only in those batches of kit from which they originate**
- **Do not detect minor deterioration of kits**
- **Need for external controls**

## External controls

- **Set of controls included from outside**
- **Put to continually monitor assay performance (validity)**
- **Positive (Borderline Reactive) & Negative**
- **Detect minor error in the assay performance**

# Why use Third Party Control/IQC/External ?

- Kit controls only calibrators, used to calculate cutoff
- Acceptance range of internal controls wide,  $PC > 0.800$
- Can detect systematic error, equipment malfunction, pipetting errors
- Internal control cannot track changes i.e. trends
  - A problem is about to occur, “near miss “ has occurred
- External control is out for range, indicative that a low positive test sample could be false negative

*Monitor changes over a period of time*

*“Smoke detector” to detect QC problem ( fire) by identifying an incorrect result*

# Why use External Control ?


- Compliance to accreditation requirements
- Basic Quality Assurance compliance
- Demonstrate competence of Staff
- Develop confidence on the lab analyzer performance

# External Controls for TTI


- Can be prepared in house: Diluted positive specimen
- Aliquots & frozen
- Available commercially: direct use for testing

# Making Suitable Dilutions

Transfer 100 ul from previous tube to next



# Selecting a Suitable Sample Dilution





# Preparation of Levy Jennings chart

- Run this weak QC sample with each assay
- Calculate E Ratio(ELISA Ratio)  

$$= \frac{\text{OD of weak QC sample}}{\text{Cut-off value}}$$
- Statistical analysis of the QC data- by plotting on LJ-graph
- Minimum 20 values required

06/02/15 *Kouale* 1/6/15 MANIPAL HOSPITAL TRANSFUSION SERVICES BANGALORE Page 1

Test Q\_HBSAG TO DETECT HEPATITIS B Lot Nr. 61084 Note 0-0.9 S/CO

Work list 01JUN15B Operator T  
 Session 50602002 Creation date 06/02/15  
 Formula: NEG+0.1 = 0.112 S/CO

Pos.	Name	OD	Mean	CV%	S/CO	Result	Note
2 A01	IQC_HBSAG	0.229			2.04	POS	
2 B01	NC	0.012					
2 C01		0.012	0.012	0.00	0.11	NEG	
2 D01	PC	3.999					
2 E01		3.999	3.999	0.00	35.71	POS	
2 F01	2821	0.011			0.10	NEG	
2 G01	2822	0.012			0.11	NEG	
2 H01	2823	0.014			0.12	NEG	
2 A02	2824	0.010			0.09	NEG	
2 B02	2825	0.013			0.12	NEG	
2 C02	2826	0.012			0.11	NEG	
2 D02	2827	0.011			0.10	NEG	
2 E02	2828	0.012			0.11	NEG	
2 F02	2829	0.012			0.11	NEG	
2 G02	2830	0.013			0.12	NEG	
2 H02	2831	0.013			0.12	NEG	
2 A03	2832	0.010			0.09	NEG	
2 B03	2833	0.013			0.12	NEG	
2 C03	2834	0.011			0.10	NEG	
2 D03	2835	0.012			0.11	NEG	
2 E03	2836	0.014			0.12	NEG	
2 F03	2837	0.013			0.12	NEG	

# Why ELISA Ratio?

---

- Why plot ELISA ratios? Why not plot OD's directly?
- Plotting E-ratio has these advantages
  - Minimizes kit to kit variation(company to company lot to lot):  
The difference in OD's and cut-off's get negated by E-ratio
  - Minimizes technician to technician variation
  - Minimizes instrument to instrument variation



# Why L-J charts?

---

- Easy to plot
- Easy to interpret
- Easy to apply Westgard rules
- Easy availability of online tools

# How to create L-J charts?

---

By using common materials

- Manually - arithmetic graph paper
- MS Excel Software in computer
- Quality software programs
- Inbuilt into Autoanalyzer



# How to create L-J chart?

By using simple statistics

- Mean




**TARGET**

- Standard Deviation (SD)

–  $\pm 1SD$

–  $\pm 2SD$

–  $\pm 3SD$



**CONTROL  
LIMITS**

- Coefficient of Variation (CV %)

# Creating L-J chart ...


- The mean and standard deviation of the control being used should be determined based on at least 20 measurements over 20 days.
- Statistical data should be determined based on >10 measurements. *(In case of perishable and less quantity control material).*

# Mean /Target Value

- **Mean** (  $\bar{x}$  ) is the sum of all the measurements ( $\Sigma$ ) divided by the number of measurements (n)

$$\text{Mean} = \bar{x} = \frac{x_1 + x_2 + x_3 + x_4 + \dots + x_N}{N}$$

N = the number of data points in the set



QC at Analytical stage

## **Standard Deviation**-Dispersion from target value

- SD quantifies the degree of dispersion of data points about the mean.

$$SD = \sqrt{\frac{(x_1 - \bar{x})^2 + (x_2 - \bar{x})^2 + (x_3 - \bar{x})^2 + \dots + (x_N - \bar{x})^2}{N}}$$



# Co-efficient of Variation %




---

- Coefficient of Variation
- Ratio of the standard deviation to the mean expressed as a percentage.

$$CV \% = \frac{\text{Standard Deviation}}{\text{Mean}} \times 100$$

- Statistical equalizer
- Performance of kits and reagents can be evaluated

# L-J chart is ready for QC monitoring



QC at Analytical stage

# Calculate Control Limits

SD is used to set limits upon which control result acceptability is determined.

$$\text{Mean} + (3 \times \text{SD}) = +3\text{SD}$$

$$\text{Mean} + (2 \times \text{SD}) = +2\text{SD}$$

$$\text{Mean} + (1 \times \text{SD}) = +1\text{SD}$$


Upper  
Control  
Limits

$$\text{Mean} - (1 \times \text{SD}) = -1\text{SD}$$

$$\text{Mean} - (2 \times \text{SD}) = -2\text{SD}$$

$$\text{Mean} - (3 \times \text{SD}) = -3\text{SD}$$

Lower  
Control  
Limits



QC at Analytical stage

# Monitoring QC Data

- Plot Levy-Jennings each day, make decision regarding acceptability
- Meet pre-established limits, usually 95% confidence interval
- Monitoring over time provides information about the ***precision*** and long-term ***accuracy*** of measurement
- Review charts at defined intervals
- Consider using Westgard Control Rules

# Monitoring QC Data

- Plot Levy-Jennings each day, make decision regarding acceptability
- Meet pre-established limits, usually 95% confidence interval
- Monitoring over time provides information about the ***precision*** and long-term ***accuracy*** of measurement
- Review charts at defined intervals
- Consider using Westgard Control Rules

# QC Results

- Random Error – factors affecting precision
  - Change in volume, temp, poor mixing, sample preparation errors, centrifugation
  - Transcription errors
  - Poor Technical skill: Lack of training, not following SOP
- Systematic error – factors affecting accuracy
  - Change in kit lot, calibrators
  - Reagent stability, degradation, deterioration
  - Inaccurate storage of reagents
  - Incorrect instrument settings

# QC Results: Systematic Variation


- Shift – six consecutive control values fall on one side of the mean, values approx equal
  - New kit lots
  - New reagents, buffers
  - Change in incubation temperature
  - New technologist, different pipetting techniques
  - Change in washer, reader

# QC Results: Systematic Variation


- Trend – six consecutive points fall in one general direction, gradual increase or decrease in value
  - Deterioration of reagents, conjugates, controls
  - Faltering equipment e,g. pipette slowly losing calibration
  - Gradual deterioration of equipment






# Shift and Trend



# WESTGARD RULES



Dr. James O Westgard



QC at Analytical stage


# Westgard Rules (James Westgard)

QC Rule	Definition	Type of error
1 - 2s	1 result $> \pm 2$ SD	Random or Systematic
1 - 3s	1 result $> \pm 3$ SD	Random or Systematic
2 - 2s	2 consecutive results $> \pm 2$ SD	Random or Systematic
4 - 1s	4 consecutive results $> \pm 1$ SD	Systematic
10x	10 consecutive results on one side of the mean	Systematic

## Westgard – 1 2s Rule

- “warning rule”
- Control result falls outside  $\pm 2SD$
- Alerts tech to possible problems
- Not cause for rejecting a run

$1_{2s}$  Rule = A warning to trigger careful inspection of the control data



# Westgard – 1 3s Rule


- If control result falls outside of  $\pm 3SD$ , rule is violated
- Run to be rejected, run regarded as Out of Control



# Westgard – 2 2s Rule

- 2 consecutive control values for the same level fall outside of  $\pm 2SD$  in the same direction
- Warning rule, requires corrective action

$2_{2s}$  Rule =  
+2SD or -2SD control limit






# Westgard – 4 1s Rule


- Four consecutive QC results for one level of control are outside  $\pm 1SD$
- Warning rule

## Control Rule Violations 4-1 S



# Westgard – 10<sub>x</sub> Rule

- Ten consecutive QC results for one level of control are on one side of the mean
- Needs attention, may not be applicable for Qualitative assay & Automation





# What to do when Control Value is out of limit?

- Rejection Rule
- “out of control”
- Stop testing
- Identify and correct problem
- Repeat testing on the samples and controls
- Do not report results until problem is solved and controls indicate proper performance


# Root Cause Analysis & CAPA

- Document all the results
- Repeat the test run when rejection occurs
- Analyze rejections & warning
- Identify the cause of error: Equipment, procedure, personnel, environment, reagents
- Resolve the problem accordingly

# Where is the Problem??



# Root Cause Analysis





# Summary

- QA in TTI: use of External Control
- Prepare LJ chart
- Apply Westgard rules for analysis
- Rejection & warnings
- Root cause analysis
- Corrective & Preventive actions
- Regular monitoring of TTI testing & continual improvement



Rotary Bangalore  &  Initiative

Thank you

ಧನ್ಯವಾದ