**TRAINING MODULE FOR
BLOOD BANK MEDICAL OFFICERS AND
LAB TECHNICIANS**

2015

***Ministry of Health and Family Welfare
Government of India***

Contributors for Gujarati translation of this manual

ISBTI Gujarat Chapter & Federation of Blood Banks of Gujarat	
Dr. Sanjay Gupta	Chairman , ISBTI Gujarat Chapter
Dr. N V Vasavada	Chairman, Federation of Blood Banks of Gujarat
Dr. Bipin Patel	Secretary, ISBTI Gujarat Chapter
Dr. Rima Kusumgar	Executive Member, ISBTI Gujarat Chapter
Dr. Nidhi Bhatnagar,	Executive Member, ISBTI Gujarat Chapter
Dr. Mamta Shah,	Department of IHBT, AHMEDABAD
Dr. Sangita Shah	Department of IHBT, AHMEDABAD
Dr. Bhargav Prajapati	Department of Transfusion Medicine, Zydus Hospital, AHMEDABAD
Dr. Amit Prajapati	Department of Transfusion Medicine, IKDRC, AHMEDABAD
Dr. Dhara Patel	Department of IHBT, AHMEDABAD
Dr. Yazdi Italia,	Zonal Chairman, ISBTI West Zone
Dr. Farzana Kothari	Executive Member, ISBTI Gujarat Chapter
Dr. Ripal Shah	Vice Chairperson, ISBTI Gujarat Chapter
Dr. Sheetal Chandak	Executive Member, ISBTI Gujarat Chapter
Mr. Chandrakant Patel	Executive Member, ISBTI Gujarat Chapter
Mr. Sandip Shah	Executive Member, ISBTI Gujarat Chapter
Mr. Govind Patel	Executive Member, ISBTI Gujarat Chapter
Mr. Navin Barot	Executive Member, ISBTI Gujarat Chapter
Mr. Jayesh Shah	Treasurer, ISBTI Gujarat Chapter

OREWORD

Shri.LovVerma, I.A.S
Secretary, Ministry of Health and Family Welfare, Govt. of India

Foreword

Though Blood Transfusion Services in India has advanced significantly over the last 50 years, the AIDS pandemic has, in fact, brought into more focus the importance of safe blood transfusions. The Blood Safety Program in India was developed as a component of the National AIDS Control Programme which has led to several advancements in terms of better policies, increased infrastructure, facilities, quality and access to blood transfusion services and promotion of voluntary blood donation.

Access to sufficient and safe blood and blood products provided within a national blood system is a vital component in achieving universal health coverage. So far, blood transfusion services were available only through established blood banks extending up to district level. In 2003, under National Health Mission these services were made available at sub-district level through Blood Storage Units at First Referral Units. These centres were aimed at meeting the requirements of blood for pregnant women requiring blood transfusions during pregnancy or labour.

Good Laboratory Practices (GLP), Good Manufacturing Practices (GMP) and quality systems are the major challenges to the organization and management of blood transfusion services. Keeping these challenges in mind and to ensure quality in blood transfusion services, 20 training institutes have been identified under NACP. These training institutes will work towards capacity building of the blood banks through standardized training curriculum for all cadres of staff in blood banks.

Dr. Sunil D. Khaparde

Deputy Director General, Basic services and Blood Safety Division, NACO

Preface

The Blood Safety Programme in India has been housed and nurtured within the National AIDS Control Organization (NACO) since the beginning of the AIDS control activities under National AIDS Control Programme (NACP). The goals and objectives of the programme are to ensure the provision of safe and quality blood even to far-flung remote areas of the country. NACO supports a network of 1137 Blood Banks in the government and charitable sectors through provision of equipment, consumables, manpower and capacity building. Nearly 70% of the country's blood requirement is met through this network. The scenario of blood banking in India owes much of its modernization due to the programmatic efforts during first three phases of the NACP.

In the fourth phase of NACP, NACO has identified specific areas for strengthening technical and service quality standards, management structures, partnership mechanisms and monitoring and evaluation systems to achieve objectives of catering to 90% of the countries blood requirements through nationally coordinated blood transfusion services.

Having recognized the need to train all cadres of staff in blood banks/ blood centres and blood storage centres, NACO/NBTC identified training centers across the country for imparting training on all aspects of blood safety since 2006. With rapid technological advances in transfusion medicine, continued training is vital for safe and quality blood supply.

This training module has been prepared with an objective of introducing uniform standards for training of blood bank staff in the subject of transfusion medicine. This module includes all aspects of blood banking for medical officers and technicians working in the blood banks. The module is designed for the better understanding and comprehension of blood bank / blood centre processes and procedures in order to improve technical and managerial skills of the trainees.

Appropriate in-service training programmes will facilitate provision of trained manpower to help the blood bank / blood centre staff to keep abreast with the latest developments, and assure quality in blood transfusion services.

Sincere appreciation for U.S. Centers for Disease Control and Prevention-Division of Global HIV/AIDS (CDC-DGHA), India Office, and Christian Medical Association of India (CMAI) for providing technical assistance and support in compilation of this document.

(Dr. Sunil D. Khaparde)

ACKNOWLEDGEMENT

The training module for Blood Bank Medical Officers and Laboratory Technicians has been developed by Blood Transfusion Services Division, NACO, Ministry of Health and Family Welfare under the guidance and active leadership of Shri. LovVerma, I.A.S., Secretary Health.

The constant encouragement of Shri K.B. Aggarwal, IAS, Joint Secretary, NACO have greatly helped in undertaking this important activity. A special thanks to Dr. ShobiniRajan, Dr. Shanoo Mishra, STI Division, NACO and other team members for their constant effort and hard work in providing direction for the guidelines.

It is commendable to note that a comprehensive document has been reviewed with the coordinated and concerted efforts of various organizations and individuals. A list of contributors is given in the document. My heartfelt thanks to all for their expertise and time towards technical review.

I extend my sincere thanks to the U.S. Centers for Disease Control and Prevention-Division of Global HIV/AIDS (CDC-DGHA), India and Christian Medical Association of India (CMAI) for providing technical assistance and support in the compilation of this document.

Dr. Sunil D. Khaparde
Deputy Director General,
Basic services and Blood Safety Division, NACO

LIST OF CONTRIBUTORS

Expert Members (In Alphabetical Order):

Dr. Atul Sonkar	Sanjay Gandhi Post Graduate Institute of Medical Science (SGPGIMS), Lucknow.
Dr. Ashish Jain	Department of Transfusion Medicine, PGIMER, Chandigarh.
Dr. Charusmita Modi	KEM Hospital, Mumbai.
Dr. Debasish Gupta	(SCTIMST), Medical College, Thiruvananthapuram
Dr. Dolly Daniel	Christian Medical College, Vellore.
Dr. Ganjendra Gupta	Santokba Duriabhji Memorial Hospital, Jaipur.
Dr. Hari Krishan Dhawan	Post Graduate of Medical Education & Research (PGIMER), Chandigarh.
Dr. Harprit Singh	NPO, (BTS), NACO, MoHFW
Dr. Jayashree Sharma	King Edward Memorial Hospital, Mumbai.
Dr. Joy Mammen	Christian Medical College, Vellore.
Dr. Kabita Chatterjee	All India Institute of Medical Science (AIIMS), New Delhi.
Dr. K. Selvarajan	Department of Transfusion Medicine, Govt. Medical College, Chennai.
Dr. Manisha Srivastava	Department of Transfusion Medicine (BMHRC), Bhopal.
Dr. M. D. Gajjar	Department of Transfusion Medicine, BJMC, Ahmedabad.
Dr. Neelam Marwaha	Post Graduate Institute of Medical Education & Research (PGIMER), Chandigarh.
Dr. Narinder Naidu	Indian Red Cross Society, Mumbai.
Dr. Pauline Harvey	U.S. Centers for Disease Control and Prevention - Division of Global HIV / AIDS.
Dr. P. Mahalingam	Christian Medical Association of India, New Delhi
Dr. R. K. Chaudhary	Sanjay Gandhi Post Graduate Institute Medical Science (SGPGIMS), Lucknow.

Dr. R. R. Sharma	Post Graduate Institute of Medical Education & Research (PGIMER), Chandigarh.
Dr. R. N. Makroo	Apollo Hospital, Mathura Road, New Delhi.
Dr. Sunder Periyavan	Transfusion Medical Centre, NIMHANS, Bangalore.
Dr. Sabitra Basu	Tata Medical Centre, Kolkata.
Dr. Sunil Rajadhyaksha	Tata Memorial Hospital, Mumbai
Dr. Sunita Upadhyaya	U.S. Centers for Diseases Control and Prevention – Division of Global HIV/ AIDS.
Dr. Sarika	Christian Medical Association of India, New Delhi.
Dr. Veena Doda	Ram Manohar Lohia Hospital, New Delhi.
Dr. Zarin S. Bharucha	Federation of Mumbai Blood Bank

Mentors

Dr. Sunil D. Khaparde,	DDG, (BTS), NACO, MoHFW
Dr. Shobini Rajan	ADG, (BTS), NACO, MoHFW
Dr. Shanoo Mishra	PO-QC, (BTS), NACO, MoHFW

ACRONYMS

ACD	Acid/Anticoagulant Citrate Dextrose
AIDS	Acquired Immuno Deficiency Syndrome
AZT/ZDV	Zidovudine
BTS	Blood Transfusion Services
CD4	Cluster Of Differentiation
CDC	Centers For Disease Control and Prevention (United States)
CMV	Cytomegalo Virus
CMAI	Christian Medical Association of India
CuSo4	Copper Sulphate
DAT	Direct Anti-Globulin Test
DHTR	Delayed Haemolytic Transfusion Reaction
DIC	Disseminated Intravascular Coagulation
EBV	Epstein Barr Virus
EDTA	EthyleneDiamineTetraacetic Acid
EIA	Enzyme Immuno Assay
ELISA	Enzyme-Linked Immuno Sorbent Assay
FFP	Fresh Frozen Plasma
Hb	Haemoglobin
HBIG	Hepatitis Immunoglobulin
HbsAg	Hepatitis B Surface Antigen
HDN	Haemolytic Disease of the Newborn
Hg	Mercury
HTC	Hospital Transfusion Committee

HTR	Haemolytic Transfusion Reaction
IAT	Indirect Antiglobulin Test
NACO	National AIDS Control Organization
NBTC	National Blood Transfusion Council
PC	Platelet Concentrate
PEP	Post Exposure Prophylaxis
PRBC	Packed Red BloodCells
PRP	Platelet Rich Plasma
PTP	Post Transfusion Purpura
QMS	Quality Management System
SAGM	Sodiumchloride, Adenine, Glucose Mannitol
SOP	Standard Operating Procedure
SBTC	State Blood Transfusion Council
TA-GVHD	Transfusion Associated Graft Versus Host Disease
TPHA	TreponemaPallidumHaemagglutination Assay
TRALI	Transfusion Related Acute Lung Injury
TTIs	Transfusion Transmitted Infections
VCTC	Voluntary Counseling and Testing Centre
VDRL	Venereal Disease Research Laboratory

CONTENTS

Sl. No	Index	Page
	List of Tables	12
	List of Figures	14
	Introduction	15
1.	Blood Donor Selection	17
1.1	Donor recruitment and retention	17
1.2	Donor Room Procedure	19
1.3	Organization of blood donation camp	40
2.	Immunohaematology (IH)	47
2.1	Basic Red Cell Serology	47
2.2	ABO and Rh blood group typing	51
2.3	ABO grouping discrepancies	57
2.4	Antiglobulin test	61
2.5	Compatibility testing	65
2.6	Quality Control of Immunohaematology reagents	70
3.	Transfusion Transmitted Infections (TTIs)	81
3.1.	Characteristics of Transfusion Transmitted Pathogens	81
3.2.	Human Immuno Deficiency Virus	83
3.3.	Hepatitis B Virus	84
3.4.	Hepatitis C Virus	85
3.5.	Syphilis	86
3.6.	Malaria	87
3.7.	Generalised Course of TTI	88
3.8.	Types of Assay	89
3.9.	Selection of Assays	93
3.10.	Blood Screening Process	93

3.11. Pre-Donation Testing	94
4. Bio safety	98
4.1. Risk to Health Care Workers (HCWs)	98
4.2. Good Laboratory Practices	100
4.3. Biomedical Waste Management	104
4.4. Disinfection	106
4.5. Autoclaving	106
5. Blood Components	109
5.1. Constituents of Blood	109
5.2. Blood Component Separation	111
5.3. Steps for Component Preparation	114
5.4. Apheresis Platelets	117
5.5. Storage and Shelf Life of Components	118
5.6. Preparation of Leucodepleted Blood Components	118
5.7. Quality Control of Blood / Blood Components	119
6. Pre transfusion issues and Bedside practices	131
6.1. Procedure for Collection of Blood Sample of the patient for pre transfusion testing	131
6.2. Procedure for arranging/requisitioning of Blood Components from the Blood Transfusion Services (BTS)	132
6.3. Administration of Blood and Blood Component	133
6.4. Disposable Equipment for administration of blood	137
6.5. Paediatric patients	137
6.6. Blood Warmers	137
6.7. Concurrent fluids and medications	138
7. Transfusion reactions	142
7.1. Classification	142
7.2. Signs and Symptoms of a transfusion reaction	143

7.3. Haemolytic Transfusion Reactions	143
7.4. Evaluation of Acute Transfusion Reaction	146
7.5. Laboratory Reactions	147
7.6. Precautions to avoid Transfusion Reactions	148
8. Quality Management of Blood Transfusion Service	151
8.1. Terms and Definitions in Quality Management System	152
8.2. Organization and Management	154
8.3. Accommodation and Environment	155
8.4. Personnel	155
8.5. Equipment	156
8.6. External services and Supply	160
8.7. Process Control	160
8.8. Identification of deviation and adverse events	161
8.9. Performance Improvement	162
8.10. Documentation and Document Control	163
8.11. Records and Test Reports	164
8.12. Internal Audit and Management Review	166
9. References	169
10. Annexures	170

LIST OF TABLES

Table 1.1	Haemoglobin Concentration and Interpretation	23
Table 1.2	Quality Control of CuSO ₄	23
Table 1.3	Conditions for Permanent Deferral	25
Table 1.4	Conditions for Temporary Deferral	26
Table 1.5	Conditions for Temporary Deferral (contd...)	27
Table 1.6	Selection of the blood bag depends on various fact	33
Table 1.7	Number of beds required as per expected no. of donors	41
Table 1.8	Donation Phase and Post Donation Phase	42
Table 2.1	Characteristics of Antibodies or Immuno-globulins	47
Table 2.2	Grading of agglutination reaction	49
Table 2.3	Reagents and Action of Potentiators	51
Table 2.4	Interpretation of forward and reverse grouping results	55
Table 2.5	Direct Antiglobulin Test (DAT)	61
Table 2.6	Indirect Antiglobulin Test (IAT)	61
Table 2.7	Infant's Blood Group and Mother's Blood Group	68
Table 2.8	Choice of Blood / Blood Components	69
Table 2.9	Platelet Concentrate (PC)	69
Table 2.10	Fresh Frozen Plasma (FFP) and Single Donor Plasma (SDP)	70
Table 2.11	Quality control of ABO reagents	73
Table 2.12	Quality control of anti-A	73
Table 2.13	Quality control of anti-B	73
Table 2.14	Quality control of anti D	74
Table 2.15	Avidity, intensity and titre of anti D	74
Table 2.16	Quality control of anti-human globulin (AHG)	75

Table 2.17	Quality control of reagents red cells	76
Table 2.18	Quality control of enzymes	76
Table 2.19	Quality control of bovine serum albumin	77
Table 3.1	Characteristics of transfusion transmitted viruses	82
Table 3.2	Window Period (in days) with various methods of detection	89
Table 3.3	Comparison of ELISA with Rapid Assay	92
Table 4.1	Modes of exposure to blood borne pathogens in the laboratory	99
Table 4.2	PEP for HBV	102
Table 4.3	HIV Post-exposure prophylaxis evaluation	103
Table 4.4	Post-exposure chemoprophylaxis	103
Table 4.5	Categories of Bio-Medical Waste	104
Table 4.6	Colour coding and type of Container for Disposal of Bio-Medical Wastes	106
Table 5.1	Components of Storage Temperature and Shelf Life	118
Table 5.2	Quality control of whole blood	119
Table 5.3	Q C of Red cell concentrate (prepared from 450 ml blood)	120
Table 5.4	Quality control of Red cell concentrate in preservative solution (Adsol/SAGM)	120
Table 5.5	Quality Control of Fresh Frozen Plasma (FFP)	121
Table 5.6	Quality control of Cryoprecipitate	122
Table 5.7	Quality control of Platelet concentrates	124
Table 5.8	Quality Control of Platelet concentrate (prepared from buffy coat)	124
Table 5.9	Quality control of Platelet Concentrate (prepared by apheresis)	124
Table 8.1	Schedules and procedures for equipment maintenance and calibration	159

LIST OF FIGURES

Figure 1	Copper Sulphate Solution	24
Figure 2	Copper Sulphate solution- Specific gravity measurement	24
Figure 3	Single Blood Bag	32
Figure 4	Double Blood Bag	32
Figure 5	Triple Blood Bag System	32
Figure 6	Quadruple Blood Bag System	33
Figure 7	Penta Blood Bag System	33
Figure 8	Cleaning the phlebotomy site with Povidone Iodine	36
Figure 9	Cleaning site with Ethyl Alcohol	36
Figure 10	Blood Sample for Blood Grouping and Cross	53
Figure 11	K3 EDTA Tube	53
Figure 12	K2 EDTA Tube	54
Figure 13	Preparation of Serial dilutions	72
Figure 14	Generalized course of TTI	88
Figure 15	Colour Coded Bins	105
Figure 16	Red Bin with Infected Waste Discard	105
Figure 17	Block Bin with General Waste Discard	105
Figure 18	Constituentsof blood products	109
Figure 19	Principle of Blood Component Separation	112
Figure 20	Whole Blood Showing Colour coded label and segments	114
Figure 21	Whole Blood with tube inserted into side slits to make it tidy during storage	114
Figure 22	Red Cell Concentrate Colour coded label and segments	114
Figure 23	Fresh Frozen Plasma Unit	114
Figure 24	Preparation of Packed Red Cells and Plasma	116
Figure 25	Preparation of Packed Red Cells in Additive Solution	116
Figure 26	Concept of Quality Management	153
Figure 27	Deviation and Adverse Events	162
Figure 28	Documentation Hierarchy	164

INTRODUCTION

Blood transfusion is an essential and life-saving support within the health care system but the safety and availability of blood and transfusion is not assured. The threats include:

- An insufficient number of voluntary non-remunerated repeat regular blood donors to ensure an adequate supply of safe blood and blood components.
- Risk of the transmission of infections such as HIV, hepatitis B and C, and syphilis through unsafe transfusion.
- Unnecessary transfusions, which needlessly expose patients to the risk of acute or delayed reactions and transfusion-transmitted infections.
- Errors in the transfusion of blood and blood components.

In order to ensure sustainability and overall safety of the entire blood transfusion process, an effective quality management system needs to be implemented.

There are many definitions of quality, one of the simplest and most appropriate being “fit for purpose”. In the context of blood transfusion, this means the setting and meeting of basic national quality standards and then continual quality improvement to ensure the safety of the transfusion process.

This Training Program is designed to improve the safety and quality of blood transfusion by capacity building in all aspects of blood transfusion.

The fundamental aim of the course is therefore to instill an understanding in the participants that for blood safety, a culture of quality has to be created in BTSs and that all staff needs to be fully aware of this culture and play their part in putting it into practice. An important part of the training course, and a measure of its effectiveness, is to ensure that course participants are aware of their central role in this process and their interest and enthusiasm for quality is aroused. One of the fundamental objectives of the courses is to help participants develop a realistic plan of action for establishing an effective quality system in their own BTS.

Blood Safety

Blood safety depends on the following factors;

- The incidence and prevalence of transfusion-transmissible infections in the blood donor population.
- The effectiveness of the blood donor education and recruitment programme and procedures for donor selection and screening, including the deferral or exclusion of unsuitable donors.
- The quality of screening of all donated blood for transfusion- transmissible infections.
- The quality of blood grouping, compatibility testing, component preparation and the storage and transportation of blood products.
- The extent to which blood and blood products are prescribed only when there is no alternative to transfusion for the particular patient.
- The reliability of the system for ensuring that patients receive blood that is compatible with their blood group, red cell antibodies and other special requirements.

It is therefore important that the quality and safety of all blood and blood products be assured throughout the process from the selection of blood donors to the administration of the product to the patient.

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Table 1.1 Hemoglobin Concentration and Interpretation

Results	Hemoglobin concentration	Interpretation
Blood drop floats	Hemoglobin < 12.5 gm/dL	Defer
Blood drop sinks	Hemoglobin ≥ 12.5 gm/dL	Accept

Note:- Drop to be observed in 15 seconds

Table 1.2 : Quality Control of CuSO₄

Date of preparation of working _____ Technician _____

Date	Specific gravity CuSO ₄	Hemoglobin estimation			Signature of Doctor
		CuSO ₄ method (Hb < or ≥ 12.5 gm/dL)	By Cell Counter	Any other Hemoglobin meter/ Color strip	

CuSO₄ (stock Solution)	
Date of preparation :	_____
Volume of stock solution prepared :	_____
Weight of dry crystal of Cuso ₄ used:	_____
Specific gravity :	_____
Technician name:	_____

Figure 1:
Copper Sulphate Solution


Figure 2: Copper Sulphate solution -
Specific gravity measurement

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Table 1.3 : Deferment of blood donation

CONDITIONS (1)	PERIOD DEFERMENT (2)
(a) Abortions	6 months
(b) History of Blood transfusion	6 months
(c) Surgery	12 months
(d) Typhoid	12 months
(e) History of malaria and duly treated	3 months (endemic) 3 years (non – endemic area)
(f) Tattoo	6 months
(g) Breast feeding	12 months after delivery
(h) Immunization (Cholera Typhoid, Diptheria, Tetanus, Plague, Gammaglobulin)	15 days
(i) Rabies vaccination	1 year after vaccination
(j) History of Hepatitis in family or close contact	12 months
(k) Immunoglobulin	12 months

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• Demethylated spirit	• Oxygen cylinder with accessories
• Betadine	• Rubber gloves
• Isopropyl swab	• First aid tray
• Cotton/Gauze swabs	• Tube striper
• Pilot Tubes : Plain and EDTA	• Electronic tube sealer
• Puncture proof container	• Needle destroyer
• Artery forceps, scissors	• Blood collceting bags
• Bp Apparatus (non-mercury) / Tourniquet	• Discard Jar with 1% Sodium Hypochlorite
• Blood bag mixer /comfortable donor couch or chair	• Adhesive plaster

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Table 1.5 : Classification of Donor Adverse Reactions

LOCAL REACTIONS RELATED TONEEDLE INSERTION
Category
Vessel Injuries Haematoma Arterial puncture Thrombophlebitis
Nerve Injuries Injury of a nerve Injury of a nerve by a haematoma
Other Complications (related to needle insertion)

Tendons injury Allergic reaction (local) Infection (local)
RARE, IMPORTANT COMPLICATIONS
Related to Vessel Injury Brachial artery pseudoaneurysm Arteriovenous fistula Compartment syndrome Axillary vein thrombosis
Accidents Accidents related to Vasovagal syncope Other kinds of accidents
Cardiovascular Reactions Angina pectoris Myocardial infarct Acute neurological condition (TIA, stroke)
Related to apheresis procedures Diffuse allergic reaction Anaphylaxis Haemolysis Air Embolus
Death
Other

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Table 1.6 : Materials required during an emergency in the post donation period

i. Oral medication ➤ Analgesic tablets (Paracetamol), Calcium and Vitamin C tablets ii. Injection ➤ Inj. Adrenaline ➤ Inj. Atropine ➤ Inj. Avil ➤ Inj. Dexamethasone ➤ Inj. Hydrocortisone ➤ Inj. Phenargan ➤ Inj. Furosemide (Lasix) ➤ Inj. Dilantin ➤ Inj. Perinorm ➤ Inj. Deriphylline ➤ Inj. Soda-bi-carbonate ➤ Inj. KCL ➤ Inj. Calcium Gluconate ➤ Inj. Dopamine ➤ Inj. Rantac ➤ Inj. Diazepam ➤ Dextrose 25% 100 ml ➤ Normal Saline 500 ml ➤ Inj. Noradranaline ➤ Inj. Mefetine ➤ Inj. Diclofenac ➤ Inj. Vitamin K ➤ Inj. Tramadol ➤ Inj. Lignocaine 2% ➤ Inj. Midazolam	iii. Antisetics ➤ Betadin, Tincture Benzoine iv. Miscellaneous ➤ Laryngoscope ➤ Laryngoscope Blade ➤ Ambu Bag with Mask ➤ Reusable oxygen reservoir ➤ Torch ➤ Airway ➤ BT-Set ➤ IV-Set ➤ Syringe ➤ E.T. 7.5mm, 7.0mm ➤ Three-way stop cock ➤ Ventimask ➤ Veinflow – 18 G, 20 G ➤ Stilette ➤ Oxygen cylinder with key ➤ Artery forceps ➤ Kidney tray (Plastic) ➤ Emesis basin ➤ Scissors ➤ Adhesive Plaster ➤ Anti Histaminic Cream ➤ Heparin and Benzyl ➤ Nicotinate ointment ➤ Smelling Salt-Spirit of Ammonia ➤ Tongue Depressor ➤ Oxygen Cylinder ➤ Infusion Set ➤ Paper Bag
---	--

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Table 1.7 : Required number of beds according to the expected number of donors

No. of Donors	No. of Couches/Beds
50-100 donors	6 to 8
100-200	8 to 10
200-300	10 to 14
Above 300 donors	14 to 16

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Apparatus and equipments required for blood donation camp

1. BP apparatus (non-mercury)
2. Stethoscope
3. Blood bags (Single – 350, Double – 350/450, Triple – 350/450, and Quadruple - 350/450)
4. Donor questionnaire, Registration cards and other essential stationery
5. Weighing device for blood donors
6. Weighing device for blood bags
7. Blood collection monitor
8. Beaker, pipette and tips
9. Tourniquet
10. Bed, sheets, draw sheets, blanket/mattress, pillows, Macintosh and arm rests
11. Lancet, blades
12. Footsteps, hand grip balls
13. Portable Hb meter with curvetters / copper sulphate working solution.
14. Glass tube / Vaccum tube for pilot sample collection.
15. Big racks for test tubes.
16. Rubber bands.
17. Spirit, betadine, cotton roll.
18. Medicated adhesive tape, micropore leucoplast roll.
19. Gel packs (coolants), anti-coagulant solution.
20. Plastic waste basket and colour coded plasticbags for waste disposal
21. Donor thank you cards.
22. Emergency medical kit
23. Insulated blood bag container's (Transpot boxes) with provision for storage temperature between 2 to 10 degrees centigrade.
24. Certificates of appreciation and honour and donor badges.

25. IEC materials, Banners and post donation standee.
26. Documentation registers, other forms and formats.
27. Soap dispensers.

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
Figure No: 10 Sample layout of blood donation camp

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
Figure: 11 Internal panoramic view of blood mobile

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Table 2.2: Tube Agglutination grading chart

Scale			Macroscopically Observed Findings
0-4	0-12		
4	12		One solid agglutinate, no free red cells detected
3	10		One or two large agglutinates, a strong reaction
2	8		Medium size agglutinates, clear background
1	5		Small agglutinates, with a lot of free red cells
+/-	3		Weak granularity in the red cell suspension
0	0		No agglutinates, an even red cell suspension

Ref: (Denise M. Harmenings, Modern Blood Banking & Transfusion Practices, 5th Ed.)

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Table 2.3 : Reagents and action of potentiators

Reagents	Action
22% Albumin	Increases dielectric constant and decreases zeta potential
LISS (Low-Ionic strength Saline)	- Decreases the ionic strength of medium (0.03 M) - Increases the rate of antibody uptake and therefore decreases the incubation time.
Enzymes (e.g. papain, ficin, bromelain)	- Removes sialic acid moieties, therefore decreases zeta potential. - Remove proteins thereby increasing interfacial tension. - Causes spicule formation therefore increases potential number of contact points. Advantage : Increases reactivity to Rh, Jk, P1, Le and antibodies Disadvantage: Removes certain antigens e.g. M, N, S Fya and Fyb.
Positively charged molecules (e.g. polybrene, PEG)	- Neutralization of sialic acid moieties thereby decreasing zeta potential. - Removes water of hydration - Increases antibody uptake because it is used as low ionic strength solution.
AHG	Cross links sensitized cells to give visible agglutination.

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Blood Group	Antigens on RBCs	Antibody in serum/plasma
A	A	Anti-B

B	B	Anti-A
O	H (no A and B)	Anti-A and Anti-B
AB	A and B	Neither anti-A nor anti-B
O _h (Bombay)	No H, A or B	Anti H, anti-A and anti-B

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Table 2.4 : Interpretation of forward and reverse grouping results

Anti-A	Anti-B	Anti-AB	Lectin A1	Lectin H	Anti-D	A1 Cells	B Cells	O Cells	Auto-Control	Blood Group
3 to 4+		3 to 4+	3 to 4+	2 to 3+	3 to 4+	-	3 to 4+	-	-	A ₁ +
3 to 4+		3 to 4+	3 to 4+	2 to 3+		-	3 to 4+	-	-	A ₁ -
3 to 4+		3 to 4+	-	3 to 4+	3 to 4+	-	3 to 4+	-	-	A ₂ +
3 to 4+		3 to 4+	-	3 to 4+	3 to 4+	1 to 2+*	3 to 4+	-	-	A ₂ +*
3 to 4+		3 to 4+	-	3 to 4+		-	3 to 4+	-	-	A ₂ -
-	3 to 4+	3 to 4+	-	3 to 4+	3 to 4+	3 to 4+	-	-	-	B+
-	3 to 4+	3 to 4+	-	3 to 4+		3 to 4+	-	-	-	B -
3 to 4+	3 to 4+	3 to 4+	3 to 4+	1 to 2+	3 to 4+	-	-	-	-	A ₁ B+
3 to 4+	3 to 4+	3 to 4+	3 to 4+	1 to 2+		-	-	-	-	A ₁ B -
3 to 4+	3 to 4+	3 to 4+	-	1 to 2+	3 to 4+	-	-	-	-	A ₂ B+
3 to 4+	3 to 4+	3 to 4+	-	1 to 2+		-	-	-	-	A ₂ B -
3 to 4+	3 to 4+	3 to 4+	-	1 to 2+	3 to 4+	1 to 2+*	-	-	-	A ₂ B+*
			-	3 to 4+	3 to 4+	3 to 4+	3 to 4+	-	-	O +
			-	3 to 4+		3 to 4+	3 to 4+	-	-	O -
			-	-	3 to 4+	3 to 4+	3 to 4+	2 to 3+	-	O _h +
			-	-		3 to 4+	3 to 4+	2 to 3+	-	O _h -

* With anti-A₁ antibodies symbol for negatives

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Adsorb the three times washed test cells with polyclonal O and B sera

(4°C for 1/2 to 1 hr)

↓

Centrifuge at 3000 rpm for 5 min

↓

Discard the supernatant

↓

Wash adsorbed cells 6 times with NS and retain the last supernatant in separate test tube

Heat elution done at 56°C for 10 min with equal volume of packed RBCs and 6% BSA

Centrifuge at 3000g for 5 min

Separate the elute in separate test tube and check in parallel with the last Supernatant wash using three unspooled A, B and O cells

Agglutination with A cells and no reaction with B or O cells suggest subgroups of A

- Similarly, subgroups of B are detected using O and A polyclonal sera

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Diagnosis	Detection
Haemolytic disease of Fetus & Newborn (HDFN)	Maternal antibodies coating fetal RBCs
Haemolytic transfusion reaction (HTR)	Recipient antibody coating donor RBCs
Auto immune Haemolytic anemia (AIHA)	Auto antibodies coating self RBCs

Table 2.6 : Indirect antiglobulin Test (IAT)
(Used to detect in-vivo sensitization of RBCs)

Application	Detection
1. Compatibility testing	Recipient antibodies with donors cells
2. Antibody screening	Recipient antibodies with screening cells
3. Antibody identification	For antibody specificity
4. Antibody titration	Rh antibody titre
5. Red cell phenotype	Weak D test, K, Jk, Fy antigens

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
Figure 12: Quality control of antisera by time estimation

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Table 2.9: Quality control of ABO reagents

Parameter	Quality requirements	Frequency of control
Appereance	No turbidity, precipitate, particle or gel formation by visual inspection	Daily
Specificity	Clear cut reactions with red cells having corresponding antigens and no reaction with negative control	Daily and with each new lot

**Table 2.10 : Quality Control of Anti-A, Anti-B and Anti-AB
(Monoclonal antisera) – Frequency of control**

Parameter	Quality requirements	Frequency of control
Appereance	No turbidity, precipitate, particle or gel formation by visual inspection	Daily
Specificity	Clear cut reactions with red cells having corresponding antigens and no reaction with negative control	Daily and with each new lot
Avidity	Repisity with which it binds to antigen- Macroscopic agglutination with 50% red cell suspension using slide test – TIME (IN SECONDS)	Daily & with each new lot
Reactivity	No immue haemolysis, Rouleax formation or prozone phenomenon	Each new lot
Potency	Undiluted serum should give 3 + reaction in test tube using 3% RBC suspension at RT	Each new lot

**Table 2.11 : Quality Control of Anti-A, Anti-B and Anti-AB
(Monoclonal antisera) – Titre & Intensity**

Antisera	Type of red cell suspension (2-5%)	Avidity time (seconds)	Titre	Intensity
Anti – A	A ₁	3 – 4	1 : 256	++++
	A ₂	5 – 6	1 : 128	++ to +++
	A ₂ B	5 – 6	1 : 64	+++
	O and B	-	-	-
Anti – B	B	3 – 4	1 : 256	+++
	A ₁ B	5 – 6	1 : 128	++++ to +++++
	O and A ₁	-	-	-
Anti – AB	A ₁	3 – 4	1 : 256	++++
	B	3 – 4	1 : 256	++++
	A ₂	5 – 6	1 : 128	+++

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Table 2.12 : Quality Control of anti D

Parameter	Quality requirements	Frequency of control
Appereance	No turbidity, precipitate or gelling	Each day
Reactivity	No Haemolysis, rouleax formation or prozone phenomenon	Each new lot
Specificity	Reactions with O positive (R ₁ r) cells and no reaction with O negative (rr) cells	Each day and each new lot

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* Negative control with rr cells should be run in parallel

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Table 2.14 : Quality control of anti-human globulin (AHG)

Parameter	Quality requirements	Frequency of control
Appereance	No turbidity, precipitate or gelling	Each day
Reactivity	No Haemolysis, rouleax formation or prozone phenomenon	Each new lot
Specificity	No agglutination of unsensitized red cells	each new lot
	Agglutination of red cells sensitized with anti D having IgG	
	Agglutination of red cells sensitized with complement binding Ab	
	Agglutination of red cells coated with C3d or C3b but no agglutination with C4 coated cells	

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Table 2.15 : Quality Control of reagents red cells

Known red cells	Anti – A	Anti – B	Anti – D
A	4+	Neg	-
B	Neg	4+	-
O positive	Neg	Neg	4+
O negative	Neg	Neg	Neg

Parameter	Quality requirements	Frequency of testing
Appearance	No Haemolysis in supernatant	Daily
Reactivity	Clear cut reactions with known antisera	Daily

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Table 2.16 : Quality control of enzymes

Parameter	Quality requirements	Frequency of testing
Reactivity	No reactivity and haemolysis with inert AB serum	Each day
Potency	2+ to 3+ agglutination of red cells sensitized with weekanti D IgG type	Each lot

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Table 2.17 : Quality control of bovine serum albumin

Parameter	Quality requirements	Frequency of testing
Appearance	No precipitation or gelling	Each day

Reactivity	No agglutination of unsensitised red cells	Each new lot
	No Haemolytic activity	Each new lot
	No prozone phenomena	Each new lot
	≥ 98% Albumin by serum electrophoresis	Each new lot
Specificity	IgG anti D should give a titre of 32 to 64 with R ₁ r red cells	Each new lot

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
Figure 14: Generalized course of Transfusion Transmissible Infections (TTIs)

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Table 3.2 : TTIS with various methods of detection

Virus	Screening marker*	Assay	Comments
HIV	Anti HIV 1& 2	Immunoassay	• Essential for effective HIV screening; screening for HIV antibody is recommended as minimum standard for blood safety • Currently the most efficacious assays are combination antigen-antibody assays • Specific detection of antibodies to both HIV-1 & HIV-2 is essential.
	HIV p24 antigen	Immunoassay	• First serological marker of HIV infection.

			• A valuable target for donation screening although viral antigen is neutralized by antibody • Screening for HIV antigen only is not appropriate as levels fall as specific antibody levels rise • HIV antigen may be detected at the same time or very soon after first detection of RNA • Currently the most sensitive HIV serological assays combine detection of both antigen (p24 antigen) and antibody (anti – HIV-1 & -2). These assays are considered to be the most effective for the serological screening of donations.
	HIV RNA	Nucleic Acid Amplification technology	• First circulating marker of HIV infection but the window between the detection of HIV RNA and HIV p24 antigen may be short • Screening for HIV RNA has been implemented in a number of countries • Value of RNA screening is related to serological screening performed and incidence of infection in donors • Evaluate increased safety vs costs and logistics
Hepatitis B	HBsAg	Immunoassay	• First serological marker of HBV infection • Significant quantities of HBsAg produced and released into circulation, the majority not associated with viral nucleic acid • Essential for effective HBV screening; screening for HBsAg recommended as minimum standard for blood safety
	Anti-HBc	Immunoassay	• Used as an additional in some countries to identify resolving infections when HBsAg has declined below detectable levels but HBV DNA may still be present. • May be the only circulating marker of infection at that point • Assays may lack specificity and specific confirmation is not available • Anti-HBs levels should be determined for all anti-HBc reactive donations to identify resolved infections • Accepted in many countries that donations that are both anti-HBc reactive and have an anti-HBs level is above 100 mIU / ml are suitable for clinical use
	HBV DNA	Nucleic acid amplification technology	• First circulating marker of HBV infection, but with limited usefulness in blood screening unless individual donation testing is performed • Virus is generally low-titre and the window between the detection of HBV DNA and HBsAg is generally very

			short
Hepatitis C	Anti-HCV	Immunoassay	- Currently the most commonly used serological marker of HCV infection - Appears in response to infection, but window period from first appearance of viral RNA may be relatively long - Screening for HCV antibody recommended as minimum standard for blood safety
	HCV antigen	Immunoassay	- First serological marker of HCV infection - A valuable target for donation screening although viral antigen is neutralized by antibody - Screening for HCV antigen only is not appropriate as levels fall as specific antibody levels rise - HCV antigen may be detected at the same time or very soon after first detection of HCV RNA - Very limited assay availability - The most sensitive HCV serological assays combine detection of both antigen and antibody. These assays are considered to be the most effective for the serological screening of donations although currently only a limited number of assays are commercially available
	HCV RNA	Nucleic acid amplification	- First circulating marker of HCV infection but the window between the detection of HCV RNA and HCV antigen may be short - Screening for HCV RNA has been implemented in a number of countries, primarily for the safety of plasma for fractionation. - Value of RNA screening is related to serological screening performed and incidence of infection in donors - Evaluate increased safety versus costs & logistics
Syphilis	Antibody to Treponema pallidum	Particle agglutination (TPHA)	- First specific serological marker of syphilis infection - Essential for effective syphilis screening for specific treponemal antibody recommended as minimum standard for blood safety
	Antibody to lipoidal antigen (reagin)	RPR	- First serological marker of syphilis infection - Do not detect syphilis specific antibody - Lack sensitivity and specificity

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Table 3.3 : Comparison of ELISA with Rapid Assay

Characteristic	ELISA	Rapid Assay
Sample	Serum, Plasma	Serum, Plasma, Whole Blood
Testing	Batch	Individual
Time required	>2 hours	<30 min
Infrastructure	Specialized set up	Minimal
Subjective variation	Minimal	Inter reader variability
Quality Control	Possible	Difficult to perform
Specificity	Better than rapid	Lower than ELISA

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Sl. No.	Errors	Resolution
1.	Insufficient volume for repeating the test in case of reactive results.	Collect 3-5 ml of blood depending upon the objective of testing and strategy of testing.
2.	Haemolysis can cause false positive or false negative result in ELISA and background noise on spot test.	Repeat sample. In case it is not possible, record “sample haemolysed” in the result.
3.	Lipaemia causes pipetting error, otherwise it causes no error in the test result.	Buy widebore pipettes and pipette tips. Record “sample lipaemic” in the result.
4.	Bacterial contamination breakdowns antibody and will affect the borderline reactive result.	Refrigerate the sample for short storage (one week), store at -200 C for long storage (>one week).
5.	Frozen /thawed sample. This may affect the borderline	Avoid freeze thawing as far as possible till the test is performed.

	positive samples	
6.	Aliquoting errors	All aliquots should be labeled. A uniform method of storing samples/aliquots should be followed in the laboratory. Organize the specimen in racks and leave a tube space empty between aliquots made from different specimens.

Table 3.5 : clerical errors

Sl. No.	Errors	Resolution
1.	Logging in specimen : may be mixing in names, wrong history, wrong number etc.	Set up the log book/record properly; collect second sample from the same patient and perform testing.
2.	Result print out : the well number on worksheet and the print out not matched properly (transcription error).	Transcribe the result from the print out sheet to the worksheet carefully. Repeat the test.
3.	Lipaemia causes pipetting error, otherwise it causes no error in the test result.	Supervisory review/vigilance of results by a second person will correct the error.
4.	Reporting to the wrong person or result communicated without post test counseling and confidentiality.	Report result to the correct person after post test counseling and maintaining the confidentiality. Reporting should be as per the “National HIV Testing Policy”.

Table 3.6 : Kit Dependent errors

Sl. No.	Errors	Resolution
1.	Used after expiration date.	Use before expiration date.
2.	Mixing of reagents from different kits/lot numbers	Use reagent for the kits for which they are made. Do not mix reagents from different.
3.	Performance characteristics (sensitivity, specificity and delta values) not satisfactory.	Do not use unsatisfactory kits.
4.	Deterioration/contamination of one or	Do not use contaminated deteriorated

	more component or reagents of the kit. (faulty transport/storage).	(control values different from those mentioned in package insert) kits. Storage at the optimal temperatures and maintain cold chain during transport of kits.
5.	Intra-lot and inter-lot variation in kit performance due to faulty manufacturing	Random monitoring of performance characteristics of the kits to ensure the quality.

Table 3.7 : Technologist dependent errors

Sl. No.	Errors	Resolution
1.	Dilutional errors	Carefully calculate the volumes of all components required for the test and perform dilutions of reconstituted components accordingly to get the desired volumes.
2.	Scratching the coated antigen from the well.	Do not touch the plate well while adding the sample.
3.	Inconsistent technique for the test and controls particularly quality control samples.	Treat each sample as well as controls in exactly the same way.
4.	Mixing reagents from different lots of kits.	Use reagents for the kits for which they are prepared for use. Do not mix reagents.
5.	Mixing of samples	The worksheets should be meticulously prepared indicating which specimen goes into which well. The internal quality control samples should be randomly placed for each run.
6.	Do not follow the SOP	Each step of SOP must be diligently followed.
7.	Wrong pipetting due to carelessness.	Supervision, vigilance and training for accurate pipetting.

Table 3.8 : Equipment based errors

Sl. No.	Errors	Resolution
1.	Use of inappropriate pipette	Use appropriate volume pipette tips to deliver

	tips	accurate volumes.
2.	Equipment not maintained as per requirement.	Maintain and calibrate the equipments as per the requirement (Refer to chapter on equipment maintenance and calibration).
3.	Washer not working satisfactorily.	The washer should be placed under annual maintenance control for optimal working.
4.	Use of wrong U.V. filters	Use the U.V. filters for reading results as per the directions in package insert.
5.	Improper reader	ELISA reader must be placed under maintenance for accurate results.
6.	Refrigerator/deep freeze not working to optimum level.	Calibrate temperature to ensure that right temperatures are maintained in refrigerator/deep freeze.

Table 3.9 : environment dependent errors

Sl. No.	Errors	Resolution
1.	High and low temperature can affect the enzyme reaction and the result	Calibration of incubators to maintain the right temperatures is very important.
2.	Drying of plates	Maintain humidity so that the plates do not dry.

Table 3.10 : Non repeatable / Non – results

Sl. No.	Causes	Resolution
1.	Mislabeling of specimen	Test 2 nd sample collected from the same patient.
2.	Specimen deterioration	Parallel testing of the sample with stored sample. Repeat testing.
3.	Borderline reactors.	Follow up and repeat testing.
4.	Carelessness of the technologist	Vigilance in the laboratory. Training of personnel.
5.	SOP not changed with the change in kit.	Review SOP and change if required particularly if new type of kits are purchased.

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Table 4.1 : Modes of exposure to blood borne pathogens in the laboratory

Procedure	HCW at Risk	Source/Modes of Transmission
Collection of blood	Laboratory technical staff / phlebotomist	• Needle stick injury • Broken specimen container • Blood contamination of hand with skin lesion / breach
Transfer of Specimen	Laboratory Technician and phlebotomist support staff and transport worker	• Contaminated exterior of the container / requisition slip • Broken container • Spill / splash of specimen
Processing of specimen	Laboratory Technical Staff	• Puncture of skin • Contamination of skin / mucous membrane from contaminated work surface • Spills / splashes of specimen • Broken specimen container • Faulty techniques • Perforated gloves
Cleaning / Washing	Laboratory support staff	• Puncture of skin • Contamination of skin from contaminated glass ware

		• Spills / splashes • Contaminated work surface
Disposal of Waste specimen transport / Mailing	Laboratory support staff, transport, postal staff	• Contact with infectious waste specially sharps, broken containers • Leaking container or packaging.

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Available chlorine concentration (labeled on bottle)	0.1% for general purpose and wiping working area	1% for spillage
Sodium hypochlorite (available chlorine 4%)	2.5 ml sodium hypochlorite + 97.5 ml water	25 ml sodium hypochlorite + 75 ml water
Sodium hypochlorite (available chlorine 5%)	2.0 ml sodium hypochlorite + 98 ml water	20 ml sodium hypochlorite + 80 ml water

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Exposed person	Treatment when source is found to be		
	HBsAg - Positive	HBsAg - Negative	Source not tested or unknown
Unvaccinated	HBIG × 1* and initiate HB vaccine	Initiate HB vaccine	Initiate HB vaccine
Previously vaccinated known responder	Test for anti – HBs titre 1. If adequate, no treatment 2. If inadequate, HB vaccine booster dose	No treatment	No treatment
Previously vaccinated known non-responder	HBIG × 2 or HBIG × 1 plus 1 dose of HB vaccine	No treatment	If known high – risk source, may treat as if source were HBsAg-positive
Response unknown	Test for anti-HBs	No	Test for anti-HBs

	1. If inadequate HBIG × 1 plus HB vaccine booster dose 2. If adequate, no treatment	treatment	1. If inadequate @, HB vaccine booster dose 2. If adequate, no treatment
--	--	-----------	--

* HBIG dose 0.06 ml/kg IM

@Adequate anti-HBsAg is > 10 SRU by RIA or positive by EIA.

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Table 4.3 : HIV Post-Exposure Prophylaxis Evalution

Exposed	Status of source (see below)			HIV status
	HIV + and low Risk	HIV + and high Risk	HIV status unknown	
Mucous membrane / non-intact skin; small volume (drops)	Consider 2-drug PEP	2-drug PEP	Consider 2-drug PEP	No PEP is required if the source blood is confirmed HIV negative
Mucous membrane / non-intact skin; large volume (major blood splash)	2-drug PEP	3-drug PEP	Consider 2-drug PEP	
Percutaneous exposure: not severe solid needle, superficial	2-drug PEP	3-drug PEP	Consider 2-drug PEP	
Percutaneous exposure: severe large bore hollow needle, deep injury, visible blood in device, needle in patient artery/vein	3-drug PEP	3-drug PEP	Consider 2-drug PEP	

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Table 4.4 : Post-Exposure Chemoprophylaxis

Basic regimen : 2 drugs (NRTIs) (4 weeks therapy)	Zidovudine (AZT/ZDV) – 300 mg twice/day is used for all types of exposure + Lamivudine (3TC) – 150 mg twice day is added to increase the effectiveness of ZDV and to prevent resistance to ZDV
Expanded Regimen: 3 drugs (2 NRTIs + PI) (4 weeks therapy)	Basic Regimen (AZT/ZDV + 3TC) + Nelfinavir 750 mg three times daily or any other boosted protease inhibitor. (for higher risk categories – consult expert)

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Table 4.5 : Categories of Bio-Medical Waste

Option	Waste Category	Treatment and disposal
Category 1	Human anatomical waste	Incineration / deep burial
Category 2	Animal waste	Incineration / deep burial
Category 3	Microbiology and Biotechnology waste	Local autoclaving / microwaving / incineration
Category 4	Waste sharps	Disinfection (autoclaving /microwaving and

		mutilation / shedding) Final disposal through authorized common BMW Treatment facility or secure landfill or concrete waste sharp pit.
Category 5	Discarded Medicines and Cytotoxic drug	Incineration / destruction drug disposal in decured and fills
Category 6	Soiled waste (items contained with blood and body fluids)	Incineration / autoclaving / microwaving
Category 7	Infectious Solid waste (Tubings, catheters, intravenous sets etc.)	Disinfection by chemical treatment autoclaving / microwaving and multilation / shredding
Category 8	Chemical Waste	Disinfect chemically and discharge into drains

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Table 4.6 : Segregation and Disposal

Sl. No.	Area	Type of waste generated	Segregation	Treatment of waste Disposal
1.	Donor screening area	Swabs	Yellow bags	Incineration
		Lancets	Puncture proof Container with Hypochlorite	Landfill/Safe pit
		Micro pipette tip	Red Bag	Shredder
		Gloves	Red Bag	Shredder
2.	Blood collection room	Swabs	Yellow Bag	Incineration
		Bag Needle	Needle destroyer (puncture proof container)	Landfill/Safe pit
		Blood Bag tubing	Red Bag	Shredder
		Gloves	Red Bag	Shredder
		Wrappers (All	Black Bag	Municipal

		Labs)		garbage
3.	Component lab	Blood bag tubing	Red Bag	Shredder
		Apheresis Kits	Red Bag	Shredder
		Gloves	Red Bag	Shredder
4.	TTI Lab	Sample tubes (glass)	Hypochlorite solution	Recycle
		Sample tubes (Plastic)	Red Bag	Shredder
		Discarded plates	Red Bag	Shredder
		Washer waste	Hypochlorite solution	Drain in sewer line
		Chemical reagents, kits controls	Hypochlorite solution	Drain in sewer line
		Reactive Bags	Autoclave/ Red Bag	Shredder
		Broken Sample tubes (Glass)	Hypochlorite solution	Land fill
		Gloves	Red Bag	Shredder
5.	Serology Lab	Sample Vials (Glass)	Hypochlorite solution	Recycle / Land fill
		Blood Bag tubing	Red Bag	shredder

Figure 15: Universal Bio Hazard Symbol

Figure 16: Colour Coded Bins

Figure 17: Red Bin with Infected Waste Discard

Figure 18: Black Bin with General Waste Discard

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Table 5.2 : Quality control of whole blood

(Ref : DGHS Technical manual)

Parameter	Quality requirement	Frequency of control
Volume	350/450 ml ± 10%	1% of all units
Anticoagulants	49/63 ml	1% of all units
PCV (Hct)	30-40%	4 units / month
Serology (HIV 1 + 2, HBsAg, HCV, MP, Syphilis)	Negative	All units
Strerility	By culture	4 units / month whichever is higher

Rest investigations as whole blood.

Table 5.3 : Quality control of Red cell concentrate
(prepared from 450 ml blood) (Ref : DGHS Technical manual)

Parameter	Quality requirement	Frequency of control
Volume	280 ml $\pm$ 40 ml	1% of all units
PCV (Hct)	70% $\pm$ 5%	Periodically
Strerility	By culture	Periodically (1% of all units)

Table 5.4 : Quality control of Red cell concentrate in
preservation solution (Adsol / SAGM) (Ref : DGHS Technical manual)

Parameter	Quality requirement	Frequency of control
Volume	350 ml $\pm$ 20 ml	1% of all units
PCV (Hct)	60% $\pm$ 5%	Periodically
Strerility	By culture	Periodically (1% of all units)

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**Table 8.1 Schedules and Produceres for equipment
maintenance and calibration**

Equipment	Performance Check
Blood collection monitor	Observe weight / volume of the first blood bag filled
Copper sulfate solution	Check specific grvity for each lot number mentioned.
Emergency power supply	Check at least once monthly
Freezers, refrigerators and platelet incubators	Monitored by temperature recording device, digital recorder or by thermometers daily. Thermographs to be reviewed weekly. Alarm system to be checked monthly.
Digital temperature recorder	Compare with thermometers daily.
Refrigerated centrifuge	Observe nd record temperature and machines speed each day of use.
General lab centrifuge	Check speed by tachometer every 6 months of use
Micro hematocrit centrifuge	Standardize before initial use, after repair or adjustments and annually
pH meter	Check pH on the day of use
Resuscitation equipment	Test flow meter and pressure gauge of O ₂ tank on day of use
Laboratory water baths	Read and record temperature daily before test.
Wayer baths for thawing Fresh Frozen Plasma	Read and record temperature daily before thawing
Shakers or rotators	Calibrate oscillations per minute daily before use.

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(I) : Sample of QC Of Components

XYZ Blood Bank

Quality control of Whole Blood and packed Red Cell

S. No.	Date	Unit No.	Type of component	Date of collection	Volume (ml)	Het (%)	WBC count (x 10 ⁶ / bag)	Plasma Hb (mg/L)	Culture	Remarks
			Outdated							
			Outdated							
			Received back							
			Received back							
			Irradiated							

Corrective measures if taken:

Signature of Technician

Signature of Doctor

(II) : Sample of QC Of Components
XYZ Blood Bank
Quality control of Platelet Components

S. No.	Date	Unit No.	Type of component	Date of collection	Volume (ml)	pH	Platelet count ($\times 10^{10}$ /unit)	WBC count ($\times 10^6$ / bag)	RBC Contamination	Swirling	Culture	Remarks

Corrective measures if taken:

Signature of Technician
Signature of Doctor

(III) : Sample of QC Of Components
XYZ Blood Bank
Quality control of Plasma Components

S. No.	Date	Unit No.	Type of component	Date of collection	Volume (ml)	PT (Sec)	APTT (sec)	F VIII (IU/unit)	FN (gm/L)	remarks
			Control for PT/APTT							

Corrective measures if taken:

Signature of Technician
 Siganature of Doctor

(IV) : Sample of QC Of Components

XYZ Blood Bank

Daily Equipment Status

Name of Equipment : _____

Manufacture : _____

Source : _____

AMC/CMC : _____

Contact No. : _____

Sl. No.	Status	Problem if any	Action taken	Remarks	Signature
1.					
2.					
3.					
4.					
5.					
6.					
7.					
8.					
9.					
10.					
11.					
12.					
13.					
14.					
15.					
16.					

Corrective measures if taken:

Signature of MO
designated

Signature of

Quality Manager

(V) : Sample of QC Of Components

XYZ Blood Bank

Calibration of Micropipettes

Date	Type of pipetee	Capacity of pipetee	Weight of aspirated water (gm)	Weight of aspirated water by standard (gm)	Result	
					Pass	Fail

Corrective mesuares if taken:

**Signature of Technician
Officer**

Signature of Mediacal

(VI) : Sample of QC Of Components
XYZ Blood Bank
Blood Discard Record

S. No.	DATE	BLOOD BAGNO	WB/PRBC	PC/BC-	PC/PRP/APH-PC	FFP / PLASMA/ CRYO / CPP	DATE OF COLLECTION	TTI-ELISA	NAT REACTIVE	EXPIRY DATE	LEAKAGE	HAEMOLYSIS	RED CELL CONT	LOW VOL	OVER VOL	DISCOLORATION	CLOTS	OTHERS

Signature of Technician
 Medical Officer

Signature of

SOP FOR L – J CHART

Subject:- Preparation of Levy – Jennigs control chart	NUMBER:- VERSION:-
LOCATION:- TTI Laboratory.	PAGES:- EFFECTIVE DATE:-
FUNCTION:- 1. To detect and investigate any variation during day to day testing carried out in TTI Laboratory. 2. To look for in between run reproducibility.	REVIEW PERIOD:- As and When required
	PREPARED BY:-
	AUTHORIZED BY:-
	DATE:-
	SIGNATURE:-

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Department of Transfusion Medicine / Blood Center:

- The sample in EDTA via (preferably vacutainer) for blood grouping and the vial must be labeled
- Requisition form and sample with discrepancy are UNACCEPTABLE
- This form will not be accepted if it is not signed or any relevant section is left blank.

Patient's Name ID No Age Sex Ward

Diagnosis.....Clinician In charge

Blood group Rh (Send fresh 2 ml sample

in EDTA vial for blood grouping if the patient has received transfusions check blood group form records as its correct information is responsibility of the doctor filling

this requisition form). Indication for Transfusion:

Pre-transfusion values:

Platelet count: _____ x 10⁵ / μ L APPT: _____ sec. PT :
_____ sec.

PTI : _____ % Serum albumin : _____ g/dL

Quantity of unit(s) required : _____

☐ Whole blood _____ ☐ Packed Cells: _____ ☐ Apheresis Platelets: _____

☐ Platelets _____ ☐ FFP / Plasma: _____ ☐ Cryoprecipitate: _____

☐ Special products

Previous Transfusion Yes ☐ No ☐

H/O Pregnancy / Abortion Yes ☐ No ☐

Adverse Reaction, if any Yes ☐ No ☐

Certify that I have personally collected the blood sample after identification of patient's ID No. and name. I have explained the necessity of component transfusion and the risks associated with it to patient / relatives and have taken informed consent.

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Name of the Hospital

CONSENT FOR THE TRANSFUSION OF BLOOD / BLOOD COMPONENTS

Patient's Name Registration No.

Ward / Bed No.

Blood transfusion is a lifesaving medical procedure. Blood can be given as 'whole blood' or as components such as : Red cells, platelets, Plasma and Cryoprecipitate.

1. I / My patient have / have been informed of the transfusion options available and expected benefits of transfusion of blood and / or components.
2. I / My patient understand/s that blood / blood components in the interest of proper medical care.
3. I / My patient understand/s that blood / blood components to be administered have been prepared and tested in accordance with rules established by National Regulation. However, there is still a very small chance that an adverse reaction can occur such as: fever with or without chills and rigor, itching and hives, which are treatable. Rarely an unpredictable life threatening event can also occur.
4. I / My patient have / has been informed that despite mandatory screening for blood borne infections such as HIV, Hepatitis B, Hepatitis C, Syphilis and Malaria, the risk of acquiring these infections is not totally eliminated.
5. I / My patient have / has had the opportunity to ask questions about transfusions, alternatives to transfusion, risk of not transfusing, the procedure to be used and the relative risks and hazards involved.
6. I / My patient believe/s that I / My patient have / has been sufficiently informed to make a decision to give consent for transfusion of blood / blood components.
7. I / My patient have / have been informed and explained the above in a language that I / My patient understand/s.

AUTHORIZATION BY PATIENT

Name of the Patient	
Signature / Thumb Impression	
Name of Witness	
Signature / thumb Impression	

Date Doctor Designation :

.....

PATIENT'S ATTENDANT / NEXT TO KIN

The patient is unable to give consent because
..... And I
..... (name / relationship to the patient),
therefore consent for the patient. I acknowledge that I have had an opportunity to
discuss this procedure, as stated above, with my physician, physician designee and
hereby consent to this procedure.

If consent is not taken please mention reason for it.

Signature of Medical Officer

**(5.1) WHOLE BLOOD / BLOOD COMPONENT / BLOOD PRODUCT
FORM**

(COMPATIBILITY REPORT)

**(form no. 1 Annexure II of Guidance document for Reporting Serious Adverse
Reactions in blood Transfusion Services, National Institute of Biologicals)**

**(Name _____ of _____ Hospoital
.....)**

(To be retained in patient's file)

1.0 PATIENT DETAILS

S. No. _____ **Date** _____
Name of Pt. _____
Age / Sex _____
Reg No. _____
Blood Group _____ **Rh** _____
Hosp. _____ **Wd** _____ **Bed** _____

2.0 PRODUCT DETAILS:

**2.1 Blood /
Components**

- 1. WB**
- 2. PRBC**
- 3. LPRBC**
- 4. PC**
- 5. PRP**
- 6. FFP**
- 7. Cryo Poor Plasma**
- 8. Cryo Precipitate**
- 9. Blood Product (Name)**

**Bag No (s) Date Blood
Bank**

- 1.**
- 2.**
- 3.**
- 4.**

Signature of Medical Officer

(5.2) WHOLE BLOOD / BLOOD COMPONENT / BLOOD PRODUCT TRANSFUSION REACTION FORM

(form no. 2 Annexure III of Guidance document for Reporting Serious Adverse
Reactions in Blood Transfusion Services, National Institute of Biologicals)

(Nam of Hospital.....)

(To be sent to Department of Transfusion Medicine / Blood Bank after transfusion)

S. No. _____ Date _____

Name of Pt. _____

Reg No. _____ Group _____ Rh _____

Blood Group _____ Rh _____

Hosp. _____ Wd _____ Bed _____

Donor Units

Blood Bag No (s.)

1.

2. Transfusion at _____ completed at _____ Rate to Transfusion _____

drops per minute Actual quantity of blood transfused _____ (ml)

CLINICAL OBSERVATION:

General condition	Pre Transfusion	During Transfusion	Post Transfusion
Pulse			
Resp.			
Temp.			
B.P.			
Rigor			
Chills			
Myalgia			
Urticaria			
Other Obseravation			

Doctor / Nurse

Note: In any case of transfusion reaction, inform the blood bank staff immediately. Send blood bag, transfusion set, post – transfusion sample (EDTA)

Signature of Medical Officer / Nurse

Blood / Components

1. WB
2. PRBC
3. LPRBC
4. PC
5. PRP
6. FFP
7. Cryo Poor Plasma
8. Cryo Precipitate
9. Blood Product (Name)

.....

(5.3) PROTOCOL FOR THE INVESTIGATION ACUTE TRANSFUSION REACTION

Annexure IV of Guidance document for Reporting Serious Adverse Reactions in Blood Transfusion Services, National Institute of Biologicals

Investigation of Acute Transfusion Reactions

- **Take immediate note and inform blood bank**
- **Seek help immediately from skilled anaesthetist or emergency team**
- **Complete the transfusion reaction form and appropriately record the following.**
 - **Type of Transfusion Reaction**
 - **Time after the start of transfusion to the occurrence of reaction**
 - **Unit No. of component transfused**
 - **Volume of the component transfused**

Send following for lab investigations:

- 1. Repeat ABO & Rh (D) Grouping**
- 2. Repeat Antibody Screen and Cross Match**
- 3. Direct Antiglobulin Test**

Send EDTA and citrated blood sample and urine sample to Hematology for:

- 4. Complete Blood Count (CBC)**
- 5. Plasma Hemoglobin**
- 6. Urine Hemoglobin**
- 7. Coagulation Screen**

Send clotted Blood sample to Biochemistry Lab for:

- 8. Renal function test (urea, creatinine and electrolytes)**
- 9. Liver function tests (bilirubin, ALT and AST)**

Send Blood culture in special blood culture bottles to Microbiology Lab.

(5.4) IMPUTABILITY LEVELS

(Sub Annexure V of Guidance document for Reporting Serious Adverse Reactions in Blood Transfusion Services, National Institute of Biologicals)

Term	Assessment scale
Definite (Certain)	When there is conclusive evidence beyond reasonable doubt

	that the adverse event can be attributed to the transfusion.
Probable (Likely)	When the evidence is clearly in favour of attributing the adverse event to the transfusion or an alternate cause.
Unlikely (Doubtful)	When the evidence is clearly in favour of attributing the adverse event to causes other than the transfusion.
Excluded	When there is conclusive evidence beyond reasonable doubt that the adverse event can be attributed to causes other than the transfusion.

(5.5) TRANSFUSION REACTION WORK-UP FORM

Annexure VI of Guidance document for Reporting Serious Adverse Reactions in
Blood Transfusion Services, National Institute of Biologicals

(Name of the Hospital)

S.No.

Name.....Hospital.....

Registration No.....Ward / Bed No.

Age / Sex..... Unit In charge

Diagnosis.....

Indication for Transfusion

Clinical Status of Patient:

Respiratory system :

Renal :

CVS:

GIT:

CNS:

Liver:

H/o Previous Transfusion, Pregnancy, Transplantation:

.....

H/o Drug intake:

Any other infusion through B.T.set:

Reaction details:

Received:

Reaction form (duly filled):

Blood Bag/Bag along with transfusion set:

Post transfusion sample:

Date/time at which blood/Blood component was transfused:

Date/time at which reaction occurred:

Date/time at which sample/reaction form were sent to blood bank

.....

TRANSFUSION REACTION WORK-UP FORM (Contd.)

(Name of the Hospital)

S.No.

Blood/BloodcomponentunitNo.

AmountofBlood/BloodComponentTransfused:

Investigation:

Identification of Patient:

Rechecking of Records:

Crossmatchfile

IssueRegister

BloodGroupingRegister

VisualExaminationofBag/Transfusionset

Supernatant of Sample:

Pre Tx Sample:

PostTxSample

BagSample

Blood Group:

PreTxSample

PostTxSample

BagSample

Direct Coombs Test (DCT):

PostTxSample

PreTxSample

Repeat cross match of Blood Bag Sample with:

Major (RT)	Major (370) AHG Phase	Minor (RT)
Pre Tx Sample		
Post Tx Sample		

TRANSFUSION REACTION WORK-UP FORM (Contd.)

(Name of the Hospital)

S.No.

Evidence of Haemolysis:

Plasma Haemoglobin Haemoglobin : Pre Tx

Serum Bilirubin: Post Tx

UrineHaemoglobin:

UrineHemosiderin :

Congulation status:

PTI:

Plateletcount:

Blood Culture (Date/time at which culture was sent):

BloodBag:

Patient:

Peripheral Blood smear (Patient sample / Blood Bag Sample):

Leishmanstain

Gramstain:

Unstainedsmear :

Blood Bag details:

TypeofBloodBag :

Lot No Date of collection :

Tube No Date of Expiry

Cross match details:

Date of Cx- match :

Emergency / Routine

Name of Technical staff who cross matched the unit

Date / time of Issue :

Interval between issue and transfusion

Where was blood kept during that interval:

Was the blood warmed before transfusion, if yes, by what method:

TRANSFUSION REACTION WORK-UP FORM (Contd.)

(Name of the Hospital)

S.No.

If Blood bag has been previously Cx-matched / issued:

Date of Cx Match	Date / time of Issue	Date of Receive Back

Donor Details

Name: Age / Sex

Address:

Phone

Date of collection : Place of collection

Name of Phlebotomist / Assitant :

Type of Donor VD / RD

Any special Investigation

Inference:

.....
.....
.....
.....
.....
.....

Signature of Consultant / Senior Resident / Junior Resident / Medical Officer

(5.7) Transfusion Reaction Reporting Form (TRRF) for Blood & Blood Components & Plasma Products

Annexure VIII of Guidance document for Reporting Serious Adverse Reactions in
Blood Transfusion Services, National Institute of Biologicals

National Institute of Biological – Indian Pharmacopoeia Commission

Ministry of Health & Family Welfare Govt. of India

HAEMOVIGILANCE PROGRAMME OF INDIA, NATIONAL INSTITUTE OF
BIOLOGICALS

(National Coordinating Centre)

For reporting of Transfusion Reactions by Healthcare Professionals

(A) PATIENT INFORMATION

Patient initials* DOB / Age in years* Blood Group* Diagnosis

Hospital Code No Hospital Admission No* Sex* F M

Date* & Time of Transfusion Date* & Time of reaction Date & Time recovery

(B) TRANSFUSION PRODUCT DETAILS*

Components	Select Components	Unit Number (Transfused)	Expiry Date of the Blood Component	Manufacture of the Blood Bag	Batch Number	Indications	1 st time / Repeat Transfusion (No. of Repeats)
Whole Blood							
Red Blood cells							
Platelets Apheresis							
Platelets pooled / RDP							
Solvent detergent (SD) Plasma							
FFP							
Cryoprecipitate							
Any other							
Plasma Products (Please Specify)	Manufacture of the Plasma Product	Batch Number	Expiry Date of the Plasma Product				

(C) NATURE OF ADVERSE REACTIONS*

Reactions	Please Tick ()
-----------	-----------------

1	Immunological Haemolysis due to ABO incompatibility															
2	Immunological Haemolysis due to other allo – antibodies															
3	Non Immunological Haemolysis															
4	Transfusion Transmitted Bacterial Infection															
5	Anaphylaxis / Hyperesensitivity															
6	Transfusion Related Acute Lung Injury (TRALI)															
7	Transfusion Transmitted Viral Infection (HBV)															
8	Transfusion Transmitted Viral Infection (HCV)															
9	Transfusion Transmitted Viral Infection (HIV - 1/2)															
10	Transfusion Transmitted Viral Infection, other (Specify)															
11	PostTransfusion Purpura															
12	Transfusion Transmitted Parasitic Infection – Malaria															
13	Transfusion Transmitted Parasitic Infection (Others specify)															
14	Transfusion Associted Graft Versus Host Disease (TAGvHD)															
15	Ferbrile Non Haemolytic Reactions (FNHTR)															
16	Transfusion Associted Dyspnea (TAD)															
17	Transfusion Associted Circulatory Overload (TACO)															
18	Other Reaction (s)															
(D)	OUTCOME OF THE ADVERSE REACTION (S)*															
	Death following the adverse reactions															
	Recoverd															
	Recovered with sequelae															
	Permanently disabled															
	Unknown															
Any other Inforamtion																
		<table border="1"> <tr> <th colspan="2">REPORTER</th> </tr> <tr> <td colspan="2">Name and professional</td> </tr> <tr> <td colspan="2">Address: _____</td> </tr> <tr> <td colspan="2">Pincode : _____</td> </tr> <tr> <td colspan="2">Email _____</td> </tr> <tr> <td colspan="2">Tel.No. (with STD code) ____</td> </tr> </table> <table border="1"> <tr> <td>CASUALTY ASSESSMENT*</td> <td>Date of this report (DD/MM/YYYY)</td> </tr> </table>	REPORTER		Name and professional		Address: _____		Pincode : _____		Email _____		Tel.No. (with STD code) ____		CASUALTY ASSESSMENT*	Date of this report (DD/MM/YYYY)
REPORTER																
Name and professional																
Address: _____																
Pincode : _____																
Email _____																
Tel.No. (with STD code) ____																
CASUALTY ASSESSMENT*	Date of this report (DD/MM/YYYY)															

(5.8) TRANSFUSION REACTION REPORTING FORM (TRRF)

Annexure VIII of Guidance document for Reporting Serious Adverse Reactions in
Blood Transfusion Services, National Institute of Biologicals

ADVICE ABOUT REPORTING

- Report adverse experience with Blood Transfusion or Blood Products administration
- Who can report?
 - Any health care professional (Doctors including Dentists, Nurses and Pharmacists)
- Where to Report?
 - Please send the completed form to the nearest Medical Colleges / Hospitals / Blood Banks / Institute under Haemovigilance Programme of India or to National Coordinating Center-Haemovigilance Programme of India, NIB
 - A list of nationwide Medical Colleges / Hospitals / Blood Banks / Institutes under Haemovigilance Programme of India is available at : <http://www.nib.gov.in>
- What happens to the submitted information
 - The causality assessment is carried out at Medical colleges / Hospitals / Blood Banks / Institutes under Haemovigilance Programme of India.
 - The information collected in Transfusion Reaction Reporting Form (TRRF) will be forwarded to National Coordinating Centre – Haemovigilance Programme of India NIB, through a software (Haemo-Vigil) developed in house by NIB's IT division. This data will be collated & analysed to identify trends, recommend best practices and interventions required to improve patient care & safety.
- These recommendations will be forwarded to IPC National Coordinating Centre, PvPI for onward transmission to Drugs Controller General (India), Central Drugs Standard Control Organization.
- These recommendations will be used to formulate safety related regulatory decisions on Blood & Blood Products Transfusion which will be communicated to various stakeholders.
- The information will be submitted to the advisory Committee of Haemovigilance Programme of India constituted by the Ministry of Health and family welfare. The Committee is entrusted with the responsibility to review any interventions that may be required.

TRANSFUSION REACTION REPORTING FORM (TRRF)

For VOLUNTARY reporting of Transfusion Reaction by

health care professionals

National Institute of Biologicals

National Coordinating Centre-Haemovigilance Programme of India

Directorate General of Health Services

Ministry of Health & Family Welfare, Government of India

A-32, Sector-62, NOIDA

<http://www.nib.gov.in>

Confidentiality


The patient's identity is held in strict confidence and protected to the fullest extent. Programme staff is not expected to and will not disclose the reporter's identity in response to a request from the public. Submission of a report does not constitute an admission that medical personnel or manufacturer or the product caused or contributed to the reactions.

TRRF can be downloaded from the websites:

www.nib.gov.in. www.ipc.gov.in. www.cdsc.nic.in

**(5.9) FLOW CHART FOR REPORTING SERIOUS ADVERSE REACTIONS
ASSOCIATED WITH BLOOD TRANSFUSION**

**Annexure IX of Guidance document for Reporting Serious Adverse Reactions in
Blood Transfusion Services, National Institute of Biologicals**

**Hemo – Vigil Software: Enter the information as per the Transfusion Reaction
Reporting form for Blood & Products for onwards transmission of data to NCC,
NIB.**

**(5.10) ISBT TABLE OF REPORTABLE SERIOUS ADVERSE
REACTIONS (SARS)**

Serious Adverse Reactions	Clinical Features	Laboratory features
Immunological: Haemolysis due to ABO incompatibility	Fever, Chills/rigors, facial flushing, chest pain, abdominal pain, back/flank pain, nausea/vomiting, diarrhea, hypotension, pallor, jaundice, oligoanuria, diffuse bleeding, dark urine, decreased haemoglobin levels. Reactions may occur within 24 hours (acute) or may not manifest for up to 28 days.	Haemoglobinuria, decreased serum haptoglobin, unconjugated hyperbilirubinaemia, increased LDH and AST levels. Blood group serology shows ABO incompatible mismatch between recipient and donor.
Immunological : Haemolysis due to other allo-antibody	As above	As above but blood group serology shows either allo-antibodies to donor red cells or auto-antibodies in the recipient
Non-immunological: haemolysis	As above	As above but due to non-immunological, possibly mechanical factors such as malfunction of a pump or blood warmer, or the use of hypotonic solutions etc.
Transfusion – transmitted bacterial infection. Note – MUST be reported	Fever, rigors and joint pain with no evidence of symptoms pre-transfusion or alternative source of infection.	Positive blood cultures from recipient and donor pack, (matching organisms) or at least one component received by the infected recipient shown to contain the agent of infection.
Transfusion – transmitted viral infection (HCV)		As above
Transfusion – transmitted viral infection HIV 1 & 2		As above
Transfusion – transmitted viral infection – other		As above
Post transfusion purpura	Bruising, severe haemorrhage, oozing wounds. Usually occurs 5 – 12 days post	Thrombocytopenia (5-12 days post transfusion) and anti-HPA antibodies present

	transfusion.	
	Fever, rash, liver, dysfunction, diarrhea. Usually occurs 1-6 weeks after transfusion. E.g. Febrile non haemolytic transfusion reactions	Pancytopenia, characteristic histological appearances on bone marrow biopsy, bone marrow hypoplasia, chimerism
Graft versus host disease	(FNHTR) where fever $\geq 39^{\circ}\text{C}$ oral or equivalent and a change of $\geq 2^{\circ}\text{C}$ from pretransfusion value, chills, rigors, headache, nausea. Usually occurs within 4 hours of transfusion and without any evidence of hemolysis, bacterial contamination or underlying condition. E.g. Transfusion associated circulatory overload (TACO) acute respiratory distress, tachycardia, increased blood pressure, acute or worsening pulmonary oedema on frontal chest x-ray, evidence of positive fluid balance. Usually occurs within 6 hours of completion of transfusion. E.g. Transfusion associated dyspnea (TAD) – respiratory distress occurring within 24 hours of transfusion but without the symptoms of TRALI, TACO or allergic reactions and not explained by any underlying condition.	

**INDUCTION TRAINING SCHEDULE FOR BLOOD BANK MEDICAL
OFFICERS
(NON BCSU)**

DAY 1					
Time	Session	Activity	Methodology	Duration	Resource person
9.00 am – 9.30 am		Registration	Group activity	30 minutes	Traning Coordinator
9.30 am – 10.10 am	Session – 1 Introduction	Activity 1 : Participant introduction	Ice breaker	15 minutes	Senior Doctor / MO blood bank
		Activity 2: Expectation of the participants	Individual Activity	5 minutes	
		Activity 3 : Pre traning assessment	Pre structured questionnair	20 minutes	
10.10 am – 10.45 am	Session 2 : Introduction to national blood safety program & National Blood Policy	Activity 1 : Over view of National blood safety program	Video	15 minutes	SACS Quality Manager
		Activity 2: National Blood Policy	Power point presentation	20 minutes	
10.45 am- 11.00 am	TEA BREAK				
11.00 am - 1.00 pm	Session 3 : Blood donor selection	Activity 1:Donor recruitment and retention	Power point presentation	30 minutes	Senior Counsellor, Senior Staff Nurse, Senior Doctor / MO Blood bank
		Activity 2: Donor room procedure & Adverse Donor Reaction	Power point presentation & Trigger video	60 minutes (video – 20 minutes, PPT – 40	

				minutes)	
		Activity 3: Organization of blood donation camp	Group activity	30 minutes	
1.00 pm – 1.45 pm	LUNCH				
1.45 pm – 3.30 pm	Session 4: Immunohaematology	Activity 1: Basic red cell serology	Power point presentation	30 minutes	Senior Doctor / MO blood bank
		Activity 2: ABO and Rh grouping & Typing	Power point presentation & Trigger video	60 minutes (video – 20 minutes, PPT – 40 minutes)	
		Activity 3: ABO grouping discrepancies	Power point presentation	15 minutes	
3.30 pm – 3.45 pm	TEA BREAK				
3.45 pm – 4.45 pm	Session 4 : Immunohaematology Contd	Activity 4: Lab demonstation HB estimation, ABO, RH grouping & typing	Lab Demo	60 minutes	Senior lab technician
DAY 2					
Time	Session	Activity	Methodology	Duration	Resource person
9.00 am – 9.30 am	Warming up	Recap DAY 1 session	Presentation	30 minutes	Rapporteur & Time Keeper
9.30 am – 10.00 am	Session – 1 Immunohematology (Contd.)	Activity 1 : Antiglobulin test	Power point presentation	30 minutes	Senior Doctor / MO blood bank, senior lab technicians
10.00		Activity 2:	Power point	30 minutes	

am – 10.30 am		Compatibilty test	presentation		
10.30 am 10.45 am	TEA BREAK				
10.45 am 11.45 am	Session 1 : Immunohematology (Contd.)	Activity 3: Group Activity on Copatibiity testing	Group activity	60 minutes	Participants
11.45 am 1.00 pm		Activity 4: QC in immunohematology	Power point presentation	70 minutes	Senior Doctor / MO blood bank,
1.00 pm – 1.45 pm	LUNCH				
1.45 pm – 2.30 pm	Session 2: Transfusion Transmissible Infection	Activity 1: Basics of TTI	Group Activity	45 minutes	Senior Doctor / MO blood bank
2.30 pm – 3.30 pm		Activity 2: Screening of TTI	Power point presentation & Trigger video	60 minutes (video – 20 minutes, PPT – 40 minutes)	Senior Doctor / MO blood bank, senior lab tachnicians
3.30 pm – 3.45 pm	TEA BREAK				
3.45 pm – 4.30 pm	Session 2: Transfusion Transmissible Infection (Contd.)	Activity 3: Laboratary Demonstation (Antiglobulin test, compatibility testing & TTI screening)	Lab Demo	45 minutes	Senior lab technician
DAY 3					

Time	Session	Activity	Methodology	Duration	Resource person
9.00 am – 9.30 am	Warming up	Recap DAY 2 session	Presentation	30 minutes	Rapporteur & Time Kepeer
9.30 am – 10.30 am	Session – 1 Bio safety & BMWM	Activity 1 : Bio safety	Power point presentation & Video	60 minutes (video – 20 minutes, PPT – 40 minutes)	Senior Doctor / MO blood bank
10.30 am – 10.45 am	TEA BREAK				
10.45 am - 11.00 am	Session 1 : Bio safety & BMWM	Activity 2: “Bin it in a minute” Game	Individual activity	15 minutes	Participants
11.00 am – 12 noon		Activity 3: Bio Medical waste management	Power point presentation	60 minutes	Senior Doctor / MO blood bank
12 noon – 1.00 pm	Session 2: Equipment management	Activity 1: Equipement & calibration	Power point presentation	60 minutes	Senior Doctor / MO blood bank, Senior lab technician
1.00 pm – 1.45 pm	LUNCH				
1.45 pm – 2.45 pm	Session 3: Blood Components	Activity 1: Blood Component	Trigger video & Power point presentation	60 minutes (Trigger video – 20 minutes, PPT – 40 minutes)	Senior Doctor / MO blood bank
2.45 pm – 3.30 pm		Activity 2: QC for components	Power point presentation	45 minutes	Senior Doctor / MO blood bank

3.30 pm – 3.45 pm	TEA BREAK				
3.45 pm – 4.00 pm	Session 3: Blood Components (contd)	Activity 3: storage & Transportation	Power point presentation	15 minutes	Senior Doctor / MO blood bank Senior lab technician
4.00 pm – 4.30 pm		Activity 4 : Administration of blood Components	Power point presentation	40 minutes	Senior Doctor / MO blood bank Senior lab technician
4.30 pm – 5.30 pm		Activity 5: Laboratory demonstration (blood Components)	Leb Demo	60 minutes	Senior lab technician
DAY 4					
Time	Session	Activity	Methodology	Duration	Resource person
9.00 am – 9.30 am	Warming up	Recap DAY 3 session	Presentation	30 minutes	Rapporteur & Time Keeper
9.30 am – 10.30 am	Session – 1 Pre transfusion issues & Bed side practices & Transfusion reaction	Activity 1 : Pre transfusion issues & Bed side practices	Power point presentation	30 minutes	Senior Doctor / MO blood bank
		Activity 2 : Transfusion reaction	Power point presentation	30 minutes	
10.30 am – 10.45 am	TEA BREAK				
10.45 am - 11.45 am	Session 2 : Roles & responsibilities of MO	Activity 1: Roles & responsibilities of MO	Power point presentation	60 minutes	Senior Doctor / MO blood bank
11.45	Session 3 : QMS in	Activity 1: “Simon	Group activity	15 minutes	By Participants

am – 1.00 pm	BTS	says” Game			
		Activity 2: Quality management & Haemovinilance	Power point presentation	60 minutes	Senior Doctor / MO blood bank
1.00 pm – 1.45 pm	LUNCH				
1.45 pm – 3.15 pm	Session 4: Record keeping documentation & legal aspects	Activity 1: : Record keeping documentation	Power point presentation	45 minutes	Senior Doctor / MO blood bank
		Activity 2: Legal aspects of Blood bank	Power point presentation	45 minutes	
3.15 pm – 3.30 pm	TEA BREAK				
3.30 pm – 4.00 pm	Session 5: Licensing and Regulation, Drugs and Cosmetic Act	Activity 1: Licensing and Regulation, Drugs and Cosmetic Act	Power point presentation	30 minutes	Senior Doctor / MO blood bank
4.00 pm – 4.30 pm	Session 6: Data analysis, Computer use and reporting	Activity 1 : Data analysis, Computer use and reporting	Power point presentation	30 minutes	
4.30 pm – 5.30 pm	Session 7: Post training assessment & wrap up	Activity 1: Post training assessment & feedback on quality of traning	Pre Structured questionnaire	30 minutes	Training Coordinator

LABORATORY DEMONSTRATION

DAY 1 : Hb estimation, ABO RH Grouping & Typing

DAY 2 : Antoglobulin test Compatibility test, TTI screening (ELISA/Rapid)

DAY 3 : Blood Components

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**INDUCTION TRAINING SCHEDULE FOR MEDICAL OFFICERS
(BCSU) & LAB TECHNICIANS**

DAY 1					
Time	Session	Activity	Methodology	Duration	Resource person
9.00 am	Registration	Registration		30 minutes	
9.30 am – 10.10 am	Session – 1 of Introduction participants	Activity 1 : Participant introduction	Ice Breaker	15 minutes	Senior Doctor / MO blood bank
		Activity 2: Expectation of the participants	Individual Activity	5 minutes	
		Activity 3 : Pre traning assessment	Individual Activity	20 minutes	
10.10 am – 10.45 am	Session 2 : Introduction to National blood safety program	Activity 1 : Over view of National blood safety program	Video Show	15 minutes	SACS Quality Manager
		Activity 2: National Blood Policy	Power point presentation	20 minutes	
10.45 am – 11.00 am	TEA BREAK				
11.00 am - 1.00 pm	Session 3 : Blood donor selection	Activity 1: Donor recruitment and retention	Power point presentation	30 minutes	Senior Senior Staff Nurse, Senior Counsellor, Doctor / MO Blood bank
		Activity 2: Donor room procedure & Adverse Donor Reaction	Trigger video, PPT	60 minutes (video – 20 minutes, PPT – 40 minutes)	
		Activity 3: Organization of blood donation	Group activity	30 minutes	

		camp			
1.00 pm – 1.45 pm	LUNCH				
1.45 pm – 3.30 pm	Session 4: Immunohaematology	Activity 1: Basic red cell serology	Power point presentation	30 minutes	Senior Doctor / MO Blood bank Senior Lab Technician
		Activity 2: ABO and Rh grouping & Typing	Trigger videos, PPT	60 minutes (video – 20 minutes, PPT – 40 minutes)	
		Activity 3: ABO grouping discrepancies	Power point presentation	15 minutes	
3.30 pm – 3.45 pm	TEA BREAK				
3.45 pm – 4.45 pm	Session 4 : Immunohaematology Contd	Activity 4: Lab demonstation HB estimation, ABO, RH grouping & typing	Lab Demonstartion	60 minutes	Senior Doctor / MO Blood bank Senior lab technician
DAY 2					
Time	Session	Activity	Methodology	Duration	Resource person
9.00 am – 9.30 am	Recap DAY 1 session	Presentation	Power point presentation	30 minutes	Rapporteur & Time Keeper
9.30 am – 10.00 am	Session – 1 Immunohematology (Contd...)	Activity 1 : Antiglobulin test	Power point presentation	30 minutes	Senior Doctor / MO blood bank, senior lab tatchnicians
10.00 am – 10.30 am		Activity 2: Compatibilty testing	Power point presentation	30 minutes	

10.30 am – 10.45 am	TEA BREAK				
10.45 am – 11.45 am	Session 1 : Immunohematology (Contd...)	Activity 3: Group Activity on Copatibiity testing	Group activity	60 minutes	Participants
1145 am – 1.00 pm		Activity 4: QC in immunohematology	Power point presentation	75 minutes	Senior Doctor / MO blood bank
1.00 pm – 1.45 pm	LUNCH				
1.45 pm – 2.30 pm	Session 2: Transfusion Transmitted Infection	Activity 1: Basics of TTI	Group Activity	45 minutes	Senior Doctor / MO blood bank
2.30 pm – 3.30 pm		Activity 2: Screening of TTI	Trigger video & PPT	60 minutes (video – 20 minutes, PPT – 40 minutes)	Senior Doctor / MO blood bank, senior lab tathnicians
3.30 pm – 3.45 pm	TEA BREAK				
3.45 pm – 4.30 pm	Session 2: Transfusion Transmitted Infection (Contd...)	Activity 3: Laboratory Demonstation (Antiglobulin test, compatibility testing & TTI screening)	Lab Demonstation	45 minutes	Senior lab technician
DAY 3					
Time	Session	Activity	Methodology	Duration	Resource

					person
9.00 am – 9.30 am	Warming up	Recap of DAY 2	Presentation	30 minutes	Rapporteur & Time Kepeer
9.30 am – 10.30 am	Session – 1 Bio safety & BMWM	Activity 1 : Bio safety	Power point presentation & Video Show	60 minutes (video – 20 minutes, PPT – 40 minutes)	Senior Doctor / MO blood bank
10.30 am – 10.45 am	TEA BREAK				
10.45 am - 11.00 am	Session 1 : Bio safety & BMWM	Activity 2: “Bin it in a minute” Game	Individual activity	15 minutes	Participants
11.00 am – 12 noon		Activity 3: Bio Medical waste management	Power point presentation	60 minutes	Senior Doctor / MO blood bank
12 noon – 1.00 pm	Session 2: Equipment management	Activity 1: Equipement & calibration	Power point presentation	60 minutes	Senior Doctor / MO blood bank, Senior lab technician
1.00 pm – 1.45 pm	LUNCH				
1.45 pm – 2.45 pm	Session 3: Blood Components	Activity 1: Blood Component	Trigger video & Power point presentation	60 minutes (Trigger video – 20 minutes, PPT – 40 minutes)	Senior Doctor / MO blood bank
2.45 pm – 3.30 pm		Activity 2: QC for blood components	Power point presentation	45 minutes	Senior Doctor / MO blood bank

3.30 pm – 3.45 pm	TEA BREAK				
3.45 pm – 4.45 pm	Session 4: Lab Demonstration	Activity 1: Laboratory demonstration (Equipment management & calibration)	Lab demonstration	60 minutes	Senior Doctor / MO blood bank Senior lab technician
DAY 4					
Time	Session	Activity	Methodology	Duration	Resource person
9.00 am – 9.30 am	Warm up	Recap DAY 3 session	Presentation	30 minutes	Rapporteur & Time Keeper
9.30 am – 11.15 am	Session – 1 Blood Components (contd...)	Activity 1 : Storage & Transportation	Power point presentation	15 minutes	Senior Doctor / MO blood bank
		Activity 2 : Administration of blood components	Power point presentation	30 minutes	Senior Doctor / MO blood bank
		Activity 3 : Apheresis (Plateletpheresis)	Power point presentation	60 minutes	Senior Doctor / MO blood bank, Senior lab technician
11.15 am – 11.30 am	TEA BREAK				
11.30 am – 12.30 pm	Session – 1 Blood Components (contd...)	Apheresis (Plateletpheresis)	Power point presentation	60 minutes	Senior Doctor / MO blood bank,

					Senior lab technician
12.30 pm – 1.00 pm	Session 2 : QC for TTI	Activity 1: Quality Control for TTI	Power point presentation	30 minutes	Senior Doctor / MO blood bank, Senior lab technician
1.00 pm – 1.45 pm	LUNCH				
1.45 pm – 2.15 pm	Session 3: Haeomovingilance	Activity 1: : Haeomovingilance	Power point presentation	30 minutes	Senior Doctor / MO blood bank, Senior lab technician
2.15 pm – 4.15 pm	Session 4: Laboratory demonstration	Activity 1: Laboratory demonstration (Component separation and Apheresis)	Lab Demo	120 minutes	Senior Doctor / MO blood bank
	TEA BREAK & WRAP UP				
DAY 5					
Time	Session	Activity	Methodology	Duration	Resource person
9.00 am – 9.30 am	Warming up	Recap DAY 4 session	Presentation	30 minutes	Rapporteur & Time Keeper
10.00 am – 10.30 am	Session – 1 Pre transfusion issues & Bed side practices & Transfusion reaction	Activity 1 : Pre transfusion issues	Power point presentation	30 minutes	Senior Doctor / MO blood bank, Senior lab technician
		Activity 2 : Bed	Power point	30 minutes	Senior

		side practices	presentation		Doctor / MO blood bank, Senior lab technician
10.30 am – 11.00 am		Activity 3 : Transfusion reaction	Power point presentation	30 minutes	Senior Doctor / MO blood bank, Senior lab technician
11.00 am – 11.15 am	TEA BREAK				
11.15 am - 12.15 pm	Session 2 : Roles & responsibilities of LT/Medical Officer	Activity 1: Roles & responsibilities of Lab Technician / Medical Officer	Power point presentation	60 minutes	Senior Doctor / MO blood bank, Senior lab technician
12.15 pm 12.45 pm	Session 3 : QMS in BTS	Activity 1: “Simon says” Game	Group activity	15 minutes	By Participants
12.45 pm – 1.00 pm		Activity 2: Quality management	Power point presentation	30 minutes	Senior Doctor / MO blood bank, Senior lab technician
1.00 pm – 1.45 pm	LUNCH				
1.45 pm – 3.15 pm	Session 4: Record keeping documentation & legal aspects	Activity 1: : Record keeping & documentation	Power point presentation	45 minutes	Senior Doctor / MO blood bank, Senior lab
		Activity 2: Legal aspects of Blood	Power point presentation	45 minutes	

		bank			technician
3.15 pm – 3.30 pm	TEA BREAK				
3.30 pm – 4.00 pm	Session 5: Licensing and Regulation, Drugs and Cosmetic Act	Activity 1: Licensing and Regulation, Drugs and Cosmetic Act	Power point presentation	30 minutes	Senior Doctor / MO blood bank, Senior lab technician
4.00 pm – 4.30 pm	Session 6: Data analysis, Computer use and reporting	Activity 1 : Data analysis, Computer use and reporting	Power point presentation	30 minutes	
4.30 pm – 5.30 pm	Session 7: Post training assessment & wrap up	Activity 1: Post training assessment & feedback on quality of training	Individual Activity	30 minutes	Training Organizer
		Activity 2: Valedictory & certificate distribution	Certificate Distribution	30 minutes	

LABORATORY DEMONSTRATION

DAY 1 : HB estimation, ABO RH Grouping & Typing

DAY 2 : Antoglobulin test Compatibility test, TTI screening

DAY 3 : Equipment Management & Calibration

DAY 4 : Component Separation and Apheresis

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
DVD – 1

Resource Materials

- 1. Training program agenda – Medical Officers, Lab Technician & Nurses.**
- 2. Training Module for Blood Bank Medical Officers & Lab Technicians**
- 3. Training Module Blood Bank Nurses**
- 4. Facilitators Guide for Blood Bank Lab Technicians & BSCU Medical Officers**
- 5. Facilitators Guide for Blood Bank Medical Officers (Non – BSCU)**
- 6. Facilitators Guide for Blood Bank Nurses**
- 7. Hand book on component Preparation and Apheresis for BSCU**
- 8. Training program power point presentation – Medical Officers, Lab Technicians & Nurses**

DVD – 2

Videos

- **Documentary**
 - (a) National Blood safety Program**
- **Trigger videos**
 - (a) Phlebotomy**
 - (b) ABO Grouping and Typing**
 - (c) Component Preparation**
 - (d) Screening of TTI**

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